

Phenolmager HT

USER GUIDE

Automated Quantitative Pathology Imaging System

2.1.0



Copyright and Disclaimers

2024 Akoya Biosciences, Inc. All rights reserved. All information in this manual and all software is subject to copyright and other intellectual property rights of Akoya Biosciences, Inc. and its affiliates and their vendors and licensors. This material may not be modified, reproduced, copied, transmitted or translated, in whole or in part, without the prior written permission of Akoya Biosciences, Inc. The information in this manual is subject to change without notice at the discretion of Akoya Biosciences, Inc.

CONTACT INFORMATION

Akoya Biosciences
100 Campus Drive, 6th Floor
Marlborough, MA 01752
United States
Phone: 855.896.8401
Fax: 855.404.0061
U.S. Email: CustomerCare@akoyabio.com
Outside U.S. Email: CustomerCareEMEA@akoyabio.com
<https://www.akoyabio.com/>

TRADEMARK

Akoya Biosciences and PhenolMager HT (registered) are trademarks of Akoya Biosciences, Inc. A Delaware corporation.

PATENT INFORMATION

Software and instruments covered by US patents 7,534,991; 7,655,898; 11,307,142; 11,879,843 and other related granted patents or pending applications throughout the world.

DISCLAIMER

This manual is published by Akoya Biosciences, Inc., 100 Campus Drive, 6th Floor, Marlborough, MA 01752 USA. Copyright 2024, Akoya Biosciences, Inc. All rights reserved, including but not limited to those rights to reproduce this publication or parts thereof. Reproduction of this publication or parts thereof, or the products it describes, by any means or in any form is expressly prohibited without the written permission of Akoya Biosciences.

LIMIT OF LIABILITY

Any errors or omission which may have occurred in this publication despite the utmost care taken in its production will be corrected as soon as possible, but not necessarily immediately upon detection. Akoya provides this publication “As Is” without warranty of any kind, either express or implied, including but not limited to the implied warranties of merchantability or fitness for a particular purpose. Some states or jurisdictions do not allow disclaimer of express or implied warranties in certain transactions; therefore, this statement may not apply to you. Akoya reserves the right to revise this publication and to make changes from time to time in the content hereof without obligation of Akoya to notify any person of such revision or changes. Further, Akoya may make modifications to the product described in this manual at any time without any obligation to notify any person of such modifications.

Contents

Preface	5
Proper Equipment Operation	5
Contact Us.....	5
Product Service and Customer Support Plans.....	5
CE	6
KC.....	6
Definition of Symbols.....	7
Safety Considerations	10
Chapter 1: Introduction.....	12
1.1. Intended Use	12
1.2. Principles of Operation.....	13
1.3. Theory of Imaging	14
1.4. Whole Slide Scanning.....	21
1.5. Multispectral Field Imaging.....	22
Chapter 2: Hardware Reference	24
2.1. Front View	25
2.2. Top View.....	27
2.3. Right-Side View	27
2.4. Left-Side Connectors.....	28
2.5. Slide Carrier Hotel	29
2.6. Slide Carrier	30
2.7. Barcode Reading	31
2.8. System Computer and Monitor	31
2.9. Specifications.....	32
Chapter 3: Operation.....	35
3.1. System Startup.....	35
3.2. System Shutdown.....	36
3.3. Inspecting Slides and Slide Carriers	36
3.4. Loading Slides into the Slide Carriers	36
3.5. Loading Slide Carriers into the Slide Carrier Hotel.....	37
3.6. Removing Slide Carriers from the Slide Carrier Hotel	38
3.7. Removing Slides from the Slide Carriers	39
Chapter 4: Software Operation	40
4.1. Software Overview	40
4.2. System Dashboard.....	43
4.3. Creating and Editing Protocols.....	48

4.4. Scanning Slides.....[69](#)

Chapter 5: Maintenance[81](#)

5.1. Power Cycling the Instrument.....[81](#)

5.2. Cleaning the Instrument Exterior.....[81](#)

5.3. Cleaning the Monitor.....[81](#)

5.4. Cleaning the Power and Communication Ports.....[82](#)

5.5. Cleaning the Slide Carriers.....[82](#)

5.6. Replacing the Fuses.....[82](#)

Appendix A: Akoya QPTIFF Specification[83](#)


Appendix B: Akoya Software & Spectra Compatibility[93](#)

Appendix C: PhenolMager HT EULA[96](#)

Preface

Section	PAGE
Proper Equipment Operation	5
Contact Us	5
Product Service and Customer Support Plans	5
CE	6
KC	6
Definition of Symbols	7
Safety Considerations	10

PROPER EQUIPMENT OPERATION

	<p>WARNING: To reduce the risk of electric shock, do not remove the instrument panels. No user serviceable parts are inside. Refer to qualified service personnel if help is required.</p> <p>Use this product only in the manner described in this manual. If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.</p> <p>AVERTISSEMENTS: Pour réduire le risque de choc électrique, ne pas retirer le couvercle. Ce produit ne contient aucune pièce pouvant être réparée par l'utilisateur. Au besoin, confier l'appareil à un réparateur qualifié.</p> <p>Ce produit ne doit être utilisé que comme décrit dans ce manuel. Si cet appareil est utilisé d'une manière autre que celle spécifiée par le fabricant, la protection fournie par l'appareil peut être entravée.</p>
--	---

CONTACT US

If you have a question about a product that is not answered in this manual or online help, or if you need assistance regarding this product, or for information about training for this product, please contact Akoya Customer Care:

Phone: +1 855-896-8401

Fax: +1 855-404-0061

Email: CustomerCare@akoyabio.com

Email (outside the US): CustomerCareEMEA@akoyabio.com

Website: www.akoyabio.com

Before calling, have the following information available for the technical representative:

- Product serial number.
- Software version (found by choosing About from the main Help menu).

PRODUCT SERVICE AND CUSTOMER SUPPORT PLANS

Akoya offers a full range of services to ensure success. From the original factory warranty through a comprehensive line of customer support plans, Akoya offers Field Service Engineers and in-house Specialists who are dedicated to supporting hardware, software, and application development needs.

Programs can include such useful services as:

- Preventive maintenance.
- Diagnostic servicing performed on-site by Akoya field service engineers or remotely via Technical Support.
- Validation performed on-site by Akoya field service engineers.
- Extended use of the Akoya Technical Support Center.
- Software updates.
- Parts, labor, and travel expense coverage.
- Other customized services upon request.

Guide Conventions

This guide contains the following information highlights and cross-references:

NOTE	<i>NOTE: Identifies important points, helpful hints, special circumstances, or alternative methods. (en français: remarque)</i>
-------------	---

CE

CE	<i>This device complies with all CE rules and requirements.</i>
-----------	--

NOTE ***Changes or modifications to this equipment not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.***

REMARQUE ***Tout changement ou modification apporté à cet instrument non expressément approuvé par l'entité responsable de la conformité peut annuler l'autorisation d'opérer l'appareil accordée à l'utilisateur.***

KOREAN CERTIFICATIONS (KC)








This device complies with MSIP (Ministry of Science, ICT, and Future Planning) EMC Registration requirements. This instrument is registered as a Class A instrument for business use only. Product sellers and users should notice that this equipment is not for household use.









A급 기기 (업무용 정보통신기기)



이 기기는 업무용으로 전자파적합등록을 한 기기이오니 판매자 또는 사용자는 이 점을 주의하시기 바라며, 만약 잘못판매 또는 구입하였을 때에는 가정용으로 교환하시기 바랍니다.

DEFINITION OF SYMBOLS

This section details the meaning, intent, and location of the labels (containing symbols) that appear on the Phenolmager® HT system and/or in the user documentation.

Symbol Symbole	Contents Description
	<p>DANGER: An imminently hazardous situation, which, if not avoided, will result in death or serious injury.</p> <p>DANGER: Situation présentant un danger imminent qui, s'il n'est pas éliminé, peut entraîner des blessures graves, voire la mort.</p>
	<p>WARNING: Caution. Refer to the User's documentation. (ISO 7000-0434B)</p> <p>AVERTISSEMENT: Attention. Se reporter à la documentation de l'utilisateur</p>
	<p>NOTE: A cautionary statement; an operating tip or maintenance suggestion; may result in instrument damage if not followed.</p> <p>REMARQUE: Énoncé indiquant une précaution à prendre, un conseil de fonctionnement ou une suggestion d'entretien; son non-respect peut provoquer des dommages à l'instrument.</p>
	<p>Hazardous voltage; risk of electric shock. (IEC 60417-6042)</p> <p>Tension dangereuse; risque de blessure par électrocution.</p>
	<p>Crush hazard. Risk of body parts, hair, jewelry, or clothing getting caught in a moving part. (ISO 3864)</p> <p>Danger d'écrasement. Faire attention que les parties corporelles, les cheveux, les bijoux ou les vêtements ne soient pas pris dans une pièce mobile.</p>
	<p>Risk of fire. (ISO 3864)</p> <p>Risque d'incendie.</p>
	<p>Risk of explosion. (ISO 3864)</p> <p>Risque d'explosion.</p>

Symbol Symbole	Contents Description
	<p>Lifting hazard. May result in injury. (ISO 3864)</p> <p><i>Levage dangereux. Peut entraîner des blessures.</i></p>
	<p>Protective ground symbol. (IEC 60417-5019)</p> <p><i>Symbole de terre de protection.</i></p>
	<p>Ground symbol. (IEC 60417-5017)</p> <p><i>Symbole de terre.</i></p>
	<p>Fuse. (IEC 60417-5016)</p> <p><i>Fusible.</i></p>
	<p>Alternating current. (IEC 60417-5032)</p> <p><i>Courant alternatif.</i></p>
	<p>On (power). (IEC 60417-5007)</p> <p><i>Marche (alimentation).</i></p>
	<p>Off (power). (IEC 60417-5008)</p> <p><i>Arrêt (alimentation).</i></p>
	<p>CE compliance mark.</p> <p><i>Marque de conformité CE.</i></p>

Symbol Symbole	Contents Description
	<p>Serial Number. (ISO 7000-2498)</p> <p><i>Numéro de série.</i></p>
	<p>Date of Manufacture. (ISO 7000-2497)</p> <p><i>Date de fabrication.</i></p>
	<p>WEEE symbol. (EN50419:2005).</p>
	<p>Unified Korea Certification Mark. (KC Mark)</p>
	<p>Signifies that the unit has passed safety tests for grounding, power line transience, and current leakage.</p> <p><i>Signifie que l'appareil a réussi les tests de sécurité pour la mise à la terre, le courant transitoire de ligne d'alimentation et la perte de courant.</i></p>
	<p>Input.</p> <p><i>Entrée.</i></p>
	<p>Output.</p> <p><i>Sortie.</i></p>
<p>Equipment labels are color coded:</p> <p><i>Les étiquettes de l'appareil sont codées couleur:</i></p>	<p>Yellow: Caution, risk of danger.</p> <p>Red: Stop.</p> <p>Blue: Mandatory action.</p> <p>Green: Safe condition or information.</p> <p><i>Jaune: Attention, danger potential.</i></p> <p><i>Rouge: Arrêter.</i></p> <p><i>Bleu: Intervention obligatoire.</i></p> <p><i>Vert: Condition sûre ou informations de sécurité.</i></p>

SAFETY CONSIDERATIONS


Safety information for the PhenolMager HT is included in this guide. Read and review all safety information before operating the PhenolMager HT instrument.

Required Training

Ensure that all personnel involved with the operation of the instrument have:



- Received instruction in general safety practices for laboratories.
- Received instruction in specific safety practices for the instrument.

DANGER, WARNINGS, AND CAUTION SIGNS

	<p><i>WARNING: Use this product only in the manner described in this manual. If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.</i></p>
---	--

Electrical Safety



The PhenolMager HT is powered by a 100-120VAC/200-240VAC, 50-60Hz (±10%) input power supply.

	<p><i>The wall outlet or the power cable connector on the left side of the instrument should be accessible after the system's installation. This enables trained service personnel to safely disconnect power from the system during servicing.</i></p> <p>Warning: Do not operate the system in an environment with explosive or flammable gases.</p>
	<p><i>Warning: DO NOT remove instrument covers. There are no user serviceable parts inside. The covers are intended to be removed only by qualified Akoya Biosciences service personnel; they are not intended to be removed during operation or for maintenance by users. Contact Technical Support at support@akoyabio.com.</i></p> <p><i>DO NOT operate the system if there has been a malfunction of the system door or slide loading components. Contact Technical Support at support@akoyabio.com if help is required.</i></p> <p><i>DO NOT operate the system in places where it may be splashed with liquid.</i></p>

For further electrical safety information, refer to the following sections:


Power Cord Selection

Contact Technical Support at support@akoyabio.com to order replacement power cords.

	<p><i>Use only the power supply cord set provided with the PhenolMager HT system. If the correct cord set for the location was not provided, contact Technical Support at support@akoyabio.com for a replacement. Do not use power supply cords with inadequate ratings.</i></p>
	<p><i>Warning: Use only a properly grounded power outlet when connecting the system to power.</i></p>


Fuses

Contact Technical Support at support@akoyabio.com to order replacement fuses.



	<p><i>Warning: The fuses in this instrument are only replaceable by trained Akoya Biosciences personnel.</i></p>
---	--

Cables and Adapters


Some cables and adapters supplied with the system have proprietary specifications.

	<p><i>Warning: Do not connect components supplied by Akoya Biosciences using unqualified cables or adapters. Contact Akoya Biosciences technical support at support@akoyabio.com to order replacement cables and adapters.</i></p>
---	--


Mechanical Safety

	<p><i>Warning: Instrument components may move during operation. Always keep body parts, hair, jewelry, and clothing away from the instrument during operation.</i></p>
	<p><i>Warning: Procedures which could result in injury may be performed only by operators who have been warned of the potential hazards and have received adequate training in performing the procedures in the safest possible manner.</i></p>

Weight Warning

	<p><i>Warning: LIFTING HAZARD. The Phenolmager HT instrument weighs 185 lbs. (84 kg). Do not move the Phenolmager HT instrument. Installing, servicing, and moving the Phenolmager HT instrument should be performed only by qualified Akoya Biosciences service personnel. Contact Technical Support at support@akoyabio.com if help is required.</i></p>
---	--

Bright Light

	<p><i>Warning: BRIGHT LIGHT HAZARD. The interior of the Phenolmager HT system includes a barcode reader with a Class 2 LED Light. Do not look into the bright light to avoid eye injury.</i></p> <p>CAUTION—CLASS 2 LASER RADIATION WHEN OPEN DO NOT STARE DIRECTLY INTO THE BEAM.</p>
---	---

Chapter 1: Introduction

Section	Page
1.1. Intended Use	12
1.2. Principles of Operation	13
1.3. Theory of Imaging	14
1.4. Whole Slide Scanning	21
1.5. Multispectral Field Imaging	22



Figure 1.1. The Phenolmager HT system.

1.1 INTENDED USE

This guide describes the use and functionality of the Phenolmager HT Automated Quantitative Pathology Imaging System. It includes operating instructions, functional descriptions, troubleshooting, illustrations, and other relevant information.

The Phenolmager HT is a multimodal digital pathology instrument that integrates both multispectral analysis and automated slide scanning that allows researchers to visualize, analyze, quantify and phenotype immune cells in situ in FFPE tissue sections and tissue microarrays (TMAs).

NOTE Akoya Biosciences Phenolmager HT Quantitative Pathology Imaging System is for research use only. Not for use in diagnostic procedures.

1.2 PRINCIPLES OF OPERATION

Akoya Biosciences PhenolMager HT is an automated imaging system for performing whole slide scans of tissue sections and microarrays (TMAs), and for acquiring multispectral (MSI) regions of interest. The system has been optimized to image samples stained with Akoya Biosciences PhenoCode™ Signature Panels (formerly Opal® multiplexed fluorescent immunohistochemistry reagent kits) and is also compatible with typical brightfield staining reagents.

The PhenolMager HT is configured to store and inventory up to 20 slide carriers, each holding up to 4 tissue slides for a total of 80 slides.

The PhenolMager HT has been designed to support a wide range of workflows:

- True whole slide scanning of slides at 1.0 $\mu\text{m}/\text{pix}$, 0.5 $\mu\text{m}/\text{pix}$, and 0.25 $\mu\text{m}/\text{pix}$.
- Analysis of multispectral whole-slide regions and fields including protein expression and phenotyping.
- Review and annotation of whole slides scans for multispectral field acquisition.
- Acquisition of multispectral fields.

To implement the full capabilities of the instrument and workflow, the PhenolMager HT system includes the following Akoya software:

PhenolMager HT

Operator-centric software for performing whole slide scans and acquiring multispectral fields of interest. This software runs on the workstation connected to the PhenolMager HT instrument.

Phenochart™

Whole-slide viewer and annotator of fluorescent and brightfield scans acquired by the PhenolMager HT. Phenochart allows the user to view the whole slide (zoom, pan, etc.), and make decisions (annotations) on next steps for the sample. Annotations in Phenochart are also used to record the workflow actions for each slide scan. Annotations include automated (inForm) analysis field requests; reviewer-requested multispectral field acquisition; and reviewer edits, approvals, and rejections. The annotation file is a fully auditable transaction log. Phenochart is freely distributed and can be used by multiple users who want to view or review slide scans taken by the PhenolMager HT.

inForm® Tissue Finder

Tissue Finder: Software typically used for the analysis of multiplexed BF and FL images. inForm supports features such as tissue classification and training, cell phenotyping, protein expression measurements, and data export. It can be run on the PhenolMager HT computer and other Microsoft® Windows® 10 computers. Additional inForm software seats beyond those that come with each PhenolMager HT system are available for purchase.

PhenolMager HT workflows range from simple 2-step procedures (e.g., acquire whole slide scan and review) to automated acquisition of regions of interest selected by the user. An additional option allows the user to go directly to whole-slide analysis using third-party software.

Example applications for the PhenolMager HT are described in the next section.

1.2.1 Example Applications

Some examples of PhenolMager HT applications include:

- Whole slide scanning of tissue samples and microarrays stained with PhenoCode Signature or Opal kits.
- Whole slide scanning of tissue samples stained with hematoxylin and eosin (H&E) and conventional immunohistochemistry (IHC) stains.
- Phenotypic analysis and protein expression of immune and cancer cells in the context of the tumor microenvironment.
- Field-based multispectral imaging (MSI) analysis of samples with up to 8 markers and counterstain with autofluorescence removal.

1.3 THEORY OF IMAGING

This section introduces some important concepts used by Akoya Biosciences PhenolMager HT imaging systems, including:

- Light
- Human perception of light intensity and color
- Light absorbance and reflection
- Fluorescence
- Photobleaching
- Dynamic range and saturation
- Multispectral Analysis
- Filter Sets for Conventional vs. Multispectral Imaging

1.3.1 Light

For purposes of this discussion, light refers to the part of the electromagnetic spectrum that can be seen by the human eye and the nearby ultraviolet and infrared wavelengths. While the physical description of light can be highly complex, we will focus on these wavelengths of light and how they interact with physical and biological materials. (See Figure 1.2.)




Figure 1.2. The electromagnetic spectrum.

1.3.2 Human Perception of Light Intensity and Color

Response to Illumination

The human eye is a highly adaptive light detector. It is significantly more sensitive in low light than in bright light. When light levels change, it takes some time for the eye to fully adjust. This is the reason people need to “dark adapt” in a darkened room before they are able to observe weak fluorescence through a microscope.

Humans can see both in very dark and very bright settings. Because the eye is so adaptable to various lighting conditions, humans are unable to quantify absolute levels of light. In any given situation, the eye has a limited ability to discriminate levels of illumination. US Department of Defense research indicates that most humans can only distinguish approximately 30 to 35 levels of gray, ranging from black to white.

The eye's response to illumination is not linear. It more closely approximates a logarithmic function. The result is that the human eye cannot see small proportional changes in brightness.

Contrasting the eye with conventional microscope imagery, any sensor that has 8-bit resolution can detect 256 levels of gray. As the number of bits of resolution increase, the number of gray levels also increases. A 12-bit sensor yields 4096 levels of gray. Digital electronic sensors are linear in response to light levels.

Ability to Distinguish Colors

While the eye is relatively poor at discriminating intensity, it is very good at distinguishing colors. Most individuals can discriminate thousands of colors. However, no two individuals see a given color in exactly the same way. The eye contains three different types of color sensors, also known as cone cells. While the arrangement of cone cells is generally standard from person to person, the ratio of each type of cone cell varies, as does their actual physical arrangement within the eye. These minor variations (along with the brain's interpretation of color) lead to the differences in perceived color between individuals.

1.3.3 Light Absorbance and Reflection

We perceive objects based on the way they transmit, absorb, and/or reflect light.

Absorbance and reflection work in tandem. Absorbance refers to the wavelengths of light that are ‘taken in by’ the objects. This means that an object that we perceive as red has absorbed all visible wavelengths of light except red. The red wavelengths are reflected back to the eye of the observer.

Transmittance refers to light emitting objects such as light sources and fluorescing or phosphorescing objects. An object we would perceive as red in transmission is one that transmits primarily red wavelengths, while absorbing or reflecting other wavelengths.

In brightfield light microscopy we observe light that passes through a specimen. Except for a few pigments and inclusions, biological specimens are essentially invisible. To impart contrast, we employ some absorbing dye, or specific optical arrangement. It is this need for contrast that led to the initial development of biological stains and stain protocols and subsequently to phase contrast and other optical contrast enhancing techniques.

Optical Density (OD) is used to measure the interaction of light with absorbing materials. The science of absorbing spectroscopy is based on the Beer-Lambert law. When absorbing images collected in brightfield are converted to OD images, the information contained in

each pixel is quantitative, as to the amount of absorbing material present. Akoya Biosciences' brightfield multispectral imagery is automatically converted to optical density at acquisition time, enabling quantitative analysis.

1.3.4 Fluorescence

Many biological and natural materials give off light of a particular color when exposed to light of another color. This property is a type of luminescence. There are two types of luminescence:

- Fluorescence refers to luminescence that occurs when the light is emitted rapidly after illumination (around one-millionth of a second).
- If the light emission takes longer than one-millionth of a second, the luminescence is called phosphorescence.

Materials that exhibit fluorescence have proven extremely useful as labels or indicators in many biological systems.

Fluorescence light emission is different than light absorption. Each fluorescent molecule generates light. Fluorescent light can be measured quantitatively because it does not interact with other materials. While fluorescence could be measured more accurately than absorbed light, a number of factors could complicate such measurements. For example, light scatters, it is affected by the local environment (such as pH), and the measurements can be affected by surrounding molecules.

Stokes Shift

When a specimen is excited with a particular (shorter) wavelength (such as blue), the specimen then shines in a different (longer) wavelength (such as red, orange, or yellow). The difference between the wavelength of the (shorter) exciting light and the wavelength of the (longer) emitted light is called the Stokes Shift, which is based on Stokes Law.

The wavelength or color used to excite the specimen (i.e., the 'excitation light') and the color the specimen glows (i.e., the 'emission light') depend on the dye involved. For any given fluorescence dye, there will be a range of excitation wavelengths that will excite fluorescence. This range of excitation wavelengths is known as the absorption spectrum. Each dye also emits across a range of wavelengths, known as the emission spectrum. Figure 1.3 contains an example of excitation and emission spectra, showing Stokes Shift and the overlap of the spectra.




Figure 1.3. Stokes Shift.

In addition, many biological materials are naturally fluorescent. This is known as autofluorescence. In particular, many vitamins, some hormones, and a variety of biological enzymes and structural proteins are naturally fluorescent. These materials often fluoresce strongly enough to interfere with specific fluorescence labeling studies.

1.3.5 Photobleaching

Because dyes can be damaged by intense light, reducing the emission signal ('photobleaching'), it is important to limit the time they are exposed to excitation light or to bright light during routine handling. Usually, blue or UV light is the most damaging. The PhenolMager HT uses an electronically-gated excitation source synchronized with its camera so the sample is only exposed to light while the camera is taking an image. Also, the PhenolMager HT front door is made of a translucent plastic that absorbs harsh blue and UV light.

NOTE *When using PhenolMager HT, avoid spending prolonged periods in the Protocol Exposures Editor while fixed on any one spot of the sample, since live image streams are being acquired.*

These steps enable repeated measurements with minimal effect on the sample.

1.3.6 Saturation and Dynamic Range

Fluorescent imaging involves measuring light emitted from a sample; the amount of light depends on the particular dye used, and the amount that is present. That, in turn, depends on the labeling and the markers involved.

In general, choose an exposure time that is long enough to measure the signals of interest without exceeding the maximum reading for the detector. Saturating the detector means the true signal level is unknown - all one can say is that it exceeded the upper measurement limit for those pixels. The exposure for each channel is chosen during assay development to balance the concerns of having good detection and avoiding saturation.

The PhenolMager HT can acquire images in either Standard or Extended dynamic range. In the latter case, readings up to 3x brighter than nominal full scale are processed without saturation. This is valuable because it allows for unexpectedly high expression in samples without any data loss.

Extended dynamic range images use more disk storage, however, because the full 12-bit camera signal is saved to disk rather than scaling it to an 8-bit data form as is done for Standard dynamic range scans.

1.3.7 Multispectral Analysis

The PhenolMager HT imaging system offers a unique solution to the problem of separating the signals from highly multiplexed samples. Multispectral analysis is based on the fact that all fluorescent materials produce a unique spectral emission. If a material is excited and the fluorescence emission examined over a range of wavelengths, then repeated for various combinations of excitation and emission wavelengths, the result is a "spectral fingerprint" which is different for each specific fluorescent material. For many fluorescent labels of biological interest, the emission spectra overlap, and may be further obscured by autofluorescence from the specimen. Multispectral imaging provides a method to distinguish between many overlapping emission spectra within the same area based on unique spectra, overcoming the limitations of conventional imaging. (See Figure 1.4.)

The result is that accurate estimates are obtained of each dye despite their overlapping spectra and the presence of background autofluorescence. The PhenolMager HT achieves

this through Akoya Biosciences' patented technologies for multispectral imaging of whole-slide scans or individual fields integrated with a simplified workflow. PhenolMager HT utilizes intelligent automation and patented technology for hands-off, error-free scanning and an optical configuration designed specifically for accurate signal measurements and unparalleled sample protection. The whole slide multispectral imaging capability creates the robust workflow since fields of view do not need to be selected. This eliminates selection bias and provides a whole-slide record for easy re-analysis of imagery with new understanding.




Figure 1.4. Unmixed multispectral image of human breast cancer tissue stained against CD4 (green), CD20 (red), CD8 (yellow), FoxP3 (orange), CD68 (purple), Cytokeratin (light blue), and DAPI (dark blue) using Opal reagents.

In general, multispectral analysis generates the spectral curves for the various fluorescent dyes or materials in a specimen. In addition, it generates a spectral curve for the autofluorescence that almost always is present to some degree. Using spectral analysis algorithms, the contributions of the individual fluorescence spectra are separated. The result is a set of images representing each spectrum that contributes to the final image. The spectral measurements consist of 8 carefully chosen filter bands, for multispectral whole slide scans; or a set of emission spectra taken with a liquid crystal tunable filter, for MSI fields. (See Figure 1.5.)




Figure 1.5. Removing auto-fluorescence with field-based multispectral imaging. A. Conventional fluorescence image of tonsil tissue stained against CD4 with Opal 520 dye. B. Overlay of all 35 layers of a multispectral image acquired using five filter cubes. C. Emission spectra of pure autofluorescence (black line), pure Opal 520 (green line), and the mixture observed in 'B' (gray dashed line). The spectral references are used to 'unmix' the contribution of autofluorescence and Opal 520 at each pixel. D. Image of unmixed Opal 520 signal extracted from B with >10-fold higher signal-to-background compared to A because the autofluorescence contribution has been separated into the component image in E.

1.3.8 Filter Sets for Conventional vs. Multispectral Imaging

Filters used for conventional fluorescence imaging are often designed to image 3-4 widely separated dyes. Typically, their excitation and emission bands have moderate spectral width, chosen to match bands where the dye is inherently most responsive. There is no attempt to measure the sample autofluorescence, or to quantify and remove it.

The PhenolImager HT provides a filter set that is optimized for whole slide multispectral imaging, with a total of 8 paired excitation and emission bands.

An alternative approach for field-based multispectral imaging makes use of a liquid crystal tunable filter (LCTF). In this case, a broad emission filter is used, and a tunable filter is engaged in the imaging path. The instrument takes pictures at several wavelengths within the emission band, to map out the full shape of the dye response. This enables analysis software to identify what dye/s are present, in what amounts, in each pixel, by spectral decomposition ("unmixing"). It also enables identifying, and removing, contributions from autofluorescence.

These two approaches are illustrated in Figures 1.6 and 1.7.




Figure 1.6. Narrowband excitation (solid line) and (dashed line) filter for conventional multispectral imaging.




Figure 1.7. Narrowband excitation (solid line) and emission long-pass emission filter (dashed line) for imaging.

Ideally, the excitation filter would match the excitation maximum of the fluorescence label being used, and the emission filter would include the emission maximum. In practical terms, the filter maxima may be slightly different from the ideal case, due to limitations of filter manufacturing and because for many dyes the Stokes shift is small, so the maxima are quite close to one another.

To image multiple dyes conventionally, one selects dyes that have very distinct excitation and/or emission response and selects filters that are narrow enough to predominantly transmit the signal of only 1 dye at a time. (See Figures 1.8 and 1.9.)




Figure 1.8. Four dyes imaged by conventional methods.




Figure 1.9. Seven dyes are too many to image fluorescence by conventional methods.

Figure 1.8 shows the emission spectra for 4 dyes (DAPI, FITC, Cy3, and Cy5.5) that are spectrally distinct and can be separated by conventional imaging. This works because there is little overlap between adjacent dyes. Excitation filters and dye response are omitted for clarity but are similarly separable.

This approach breaks down when more dyes are present, or when it is important to account for the effects of autofluorescence. Figure 1.9 shows the spectra of 7 dyes, and with this many dyes there is no way to isolate their signals using conventional imaging techniques, or to account for autofluorescence.

1.4 WHOLE SLIDE SCANNING

PhenolMager HT scans slides using the following process:

Color Overview

PhenolMager HT takes a low power color overview of all four slides in each carrier, including the label for each slide. This initial step is performed regardless of whether a fluorescence or brightfield protocol is used and identifies the presence of slides in the carrier and capture their labels.

Coverslip Finding

For each slide, PhenolMager HT will then find the coverslip using this overview scan. The coverslip defines the potential scan area.

For fluorescence protocols, this area can be further restricted by making a closed loop with a red, green, or blue Sharpie® marker. If closed-loop markup is present, the system will only scan within the loop. This is useful if the tissue is faintly stained, punctate, or is a highly fluorescing non-tissue material (PAP pen, for example).

For brightfield protocols, the coverslip will define the potential scan area. Closed-loop markup is not available for brightfield scans.

Slide Height Finding

PhenolMager HT engages specialized height-sensing optics to measure the top of the coverslip at up to 9 locations. This gives an initial focus estimate based on the expected coverslip thickness.

Fluorescence Overview (fluorescence protocols only)

If using a fluorescence protocol for this slide, PhenolMager HT will take another overview within the coverslip (or closed-loop markup), this time in fluorescence.

Tissue Finding

Using the corresponding overview, PhenolMager HT will automatically detect the sample on the slide. The resulting area will be scanned. If requested within the protocol, PhenolMager HT will scan the entire area within the coverslip (or closed-loop markup).

Focus Finding

PhenolMager HT will measure focus at multiple points on the tissue to determine the best focus. It uses the sample map from the previous step to choose the measurement grid and continues until the grid is fully measured.

If the measured tissue height is irregular, PhenolMager HT increases the grid density and takes more readings until it finds the readings are regular at the new, finer scale.

All focus measurements include a dust-rejection algorithm, and an overall consistency check is applied to further reduce the likelihood of dust-induced focus errors.

Scanning

PhenolMager HT then scans the slide.

Brightfield scans are conventional color scans that have been color and background corrected.

Fluorescence scans are multi-layered, with one layer for each filter chosen. To avoid photobleaching, the system uses a pulsed LED so the sample is only exposed to light during the time the camera is taking a picture.

Completed scans can be opened in Phenochart, the whole slide viewing application.

NOTE *Some dyes narrowly express in a single filter. Other times, dyes may express in multiple filters and can appear in more than one layer in a fluorescence scan. This can occur when using PhenolMager HT's broad filters like the Cy3 or Texas Red filters with Opal 570 dye, or when using non-Akoya dyes. If the sample is highly multiplexed, multiple dyes may appear in the same channel.*

Spectral Unmixing

When spectral unmixing is requested, PhenolMager HT unmixes the scan in the background, while it turns to the next slide. This produces a new scan file containing the unmixed imagery.

Unmixing is highly optimized, and typically completes while the next slide is still underway. Thus, the tasks are fully nested and overall system throughput is kept very high.

1.5 MULTISPECTRAL FIELD IMAGING

In some workflows, the PhenolMager HT is used to acquire multispectral fields of a sample, after it has completed a whole slide scan of the tissue.

1.5.1 Configuring the Multispectral Fields

Multispectral imagery is acquired using the following process:

- Selection of Multispectral Fields
- Fields for multispectral imaging are selected on a previously scanned slide. Using Phenochart, select individual fields or regions of interest. If desired, also select fields using a trained inForm algorithm. See the Phenochart User Guide for more information.

1.5.2 Processing the Multispectral Fields

After configuring the slide for field acquisition, PhenolMager HT will perform the following actions:

1. Color Overview

PhenolMager HT takes a low-power color overview of all four slides in each carrier, including the label for each slide. This initial step is performed regardless of whether a fluorescence or brightfield protocol is used and identifies the presence of slides in the carrier.

2. Slide Registration

Using the color overview along with the slide's original overview, PhenolMager HT accounts for any shift or rotation of the slide to ensure that the multispectral region locations are accurate. The slide edges are used to account for any rotation or horizontal shift. The coverslip edges are used to account for vertical shifting.

3. Slide Height Finding

PhenolMager HT engages specialized height-sensing optics to measure the top of the coverslip at up to 9 locations. This gives an initial focus estimate based on the expected coverslip thickness.

4. Acquisition of Multispectral Fields


PhenolMager HT then travels to each multispectral region site, autofocuses, corrects for exposure if requested, and acquires the multispectral image.

Multispectral imagery can then be viewed and analyzed using inForm.

Chapter 2: Hardware Reference

This chapter identifies and describes the Phenolmager HT system hardware and lists the technical specifications.

Section	Page
2.1. Front View	25
2.2. Top View	27
2.3. Right-Side View	27
2.4. Left-Side Connectors	28
2.5. Slide Carrier Hotel	29
2.6. Slide Carrier	30
2.7. Barcode Reading	31
2.8. System Computer and Monitor	31
2.9. Specifications	32

	<p><i>Warning: Lifting Hazard. Do not move the Phenolmager HT instrument. Installing, servicing, and moving the Phenolmager HT instrument should be performed only by qualified Akoya Biosciences service personnel. Contact Technical Support at support@akoyabio.com if help is required.</i></p>
--	---

2.1 FRONT VIEW

The front view of the Phenolmager HT instrument with the system door open. (See Figure 2.1.) The parts identified are described in Table 2.1.




Figure 2.1. Phenolmager HT front view—system door open.

Table 2.1. Front View—hardware components

Part	Description
Status Progress Display	Each light represents an individual slide and indicates the progress of the slide currently being scanned.
Slide Carrier Hotel	Houses up to 20 slide carriers. See Slide Carrier Hotel for more detailed information.
Door Sensor	Opens or closes the System Door when a hand is placed in front of the sensor.

The front view of the Phenolmager HT instrument with the system door closed. (See Figure 2.2.) The parts identified are described in Table 2.2.




Figure 2.2. PhenolMager HT front view—system door closed.

Table 2.2. Front View (system door closed)—hardware components

Part	Description
Slide Carrier Status Lights	Each Slide Carrier Status Light represents an individual slide carrier and indicates the processing status of each Slide Carrier. See Table 2.7 for the color codes of the slide carrier status lights.
System Door	When closed, covers the Slide Carrier Hotel and slide processing can occur. When open, reveals the Slide Carrier Hotel.

	<p><i>Warning: Do not operate the system if there has been a malfunction of the system door or slide loading components. Contact Technical Support at support@akoyabio.com if help is required.</i></p>
--	---

2.2 TOP VIEW

The top view of the PhenolMager HT instrument is seen in Figure 2.3 and the Cover Handle is described in Table 2.3.




Figure 2.3. Phenolmager HT - top view.

Table 2.3. Top View—hardware components

Part	Description
Cover Handle	Used by only service personnel to remove instrument covers during service and installation.

	<p><i>Warning: DO NOT remove instrument covers. There are no user serviceable parts inside. The covers are intended to be removed by qualified service personnel only; they are not intended to be removed during operation or for maintenance by users. Contact Technical Support at support@akoyabio.com if help is required.</i></p>
--	---

2.3 RIGHT-SIDE VIEW

Figure 2.4 shows the right-side of the Phenolmager HT instrument. The System Power Switch identified is described in Table 2.4.




Figure 2.4. Phenolmager HT – right-side view.

Table 2.4. Right-Side View—hardware components

Part	Description
System Power Switch	Turns the Phenolmager HT instrument ON or OFF.

2.4 LEFT-SIDE CONNECTORS

Figure 2.5 identifies the connectors on the bottom left side of the Phenolmager HT instrument. The connectors identified are described in Table 2.5.

	<p><i>Warning: Appliance inlet is a disconnecting device. Place device or equipment in a manner so that disconnecting device is accessible at all times.</i></p> <p><i>Use only the power supply cord set provided with the Phenolmager HT system. If the correct cord set for the location was not provided, contact Akoya Biosciences Technical Support at support@akoyabio.com for a replacement.</i></p> <p><i>Use only a properly grounded power outlet when connecting the system to power.</i></p>
--	---




Figure 2.5. Phenolmager HT left-side connectors.

Table 2.5. Left-Side View—connection components

Part	Description
Power Connector	Connects to the power cord to provide power to the Phenolmager HT instrument.
USB 2.0 Connector	Connects a USB 2.0 cable to a USB 2.0 port on the System Computer Connectors.
USB 3.0 Hub Connector	Connects a USB 3.0 HUB to a USB 3.0 port on the System Computer Connectors.

2.5 SLIDE CARRIER HOTEL

The Slide Carrier Hotel houses up to 20 Slide Carriers and is visible when the System Door opens. Figure 2.6 shows a closeup of the Slide Carrier Hotel and identifies its components. The hardware and color components are described in Table 2.6 and Table 2.7 respectively.

	<p><i>Warning: Do not operate the system if there has been a malfunction of the system door or slide loading components. Contact Technical Support at support@akoyabio.com if help is required.</i></p>
---	---














Figure 2.6. Slide Carrier Hotel.

Table 2.6. Slide Carrier Hotel—hardware components

Part	Description
Slide Carrier Status Lights	Indicates the processing status of each slide carrier. The color codes associated with the slide carrier status lights are listed in Table 2.7.
Hotel Slot	Each hotel slot holds one Slide Carrier. The hotel slots are numbered 1-20, starting at the top of the hotel.
Slide Carrier	Houses up to 4 slides. See the Slide Carrier section for more details.

Table 2.7. Slide Carrier Status Lights—color codes

Color	Status
 None	Slide Carrier Hotel is empty.

Color	Status
 White	Initial state of slot after slide carrier is inserted and no rules have been applied via software.
 Solid	Hotel slot is malfunctioning.
 Blinking	Slide carrier is not completely inserted.
 Solid	Awaiting Processing.
 Blinking	Prioritized and awaiting processing.
 Solid Yellow	Processing instructions are either incomplete or invalid.
 Blinking	Processing.
 Solid Green	Processing Complete.
 Orange	Processing complete but a processing error occurred.

2.6 SLIDE CARRIER

Slide Carriers are an integral component of PhenolMager HT. The Slide Carrier holds up to 4 microscope slides and helps protect slides from damage. The microscope slide positions are labeled 1 to 4 and correspond to slide identification in the PhenolMager HT software.

Figure 2.7 shows a closeup of a Slide Carrier and identifies its components. The identified components are described in Table 2.8. Table 2.9 provides details of each slide format.



Figure 2.7. Phenolmager HT Slide Carrier.

Table 2.8. Phenoimager HT Slide Carrier—hardware components


Part	Description	
Slide Position	Each slide position is populated with one slide. Slide position 1 is on the far left; slide position 4 is on the far right.	
Slide	Each microscope slide is manually loaded into the slide positions. Table 2.9 identifies the details of the slide format.	
Insert Indicator		An icon that indicates the side of the slide carrier to be inserted into a Hotel Slot.
Tab Cover	Holds the spring-loaded tabs in place.	
Spring-Loaded Tab	When a slide is inserted into a slide slot, it is gently placed against a spring-loaded tab.	
Carrier Handle	The side of the slide carrier to hold when inserting and removing a slide carrier from a Hotel Slot. A unique number is printed on top of the handle for slide carrier identification purposes.	

Table 2.9. Details of slide formats

Part	Width	Height	Thickness
Metric	25.0 ± 1.0 mm	75.0 ± 1.0 mm	1.00 ± 0.10 mm
English	25.4 mm	76.2 mm	1.00 ± 0.10 mm

2.7 BARCODE READING

The Phenolmager HT reads barcodes on the slide labels by taking a picture of the label and then decoding it.

2.8 SYSTEM COMPUTER AND MONITOR

The Phenolmager HT system includes a widescreen monitor and a computer pre-installed with Phenolmager HT, Phenochart, and inForm software. A wireless keyboard and mouse are also included.

The Phenolmager HT computer system connectors are seen in Figure 2.8. The connectors identified are described in Table 2.10.




Figure 2.8. System Computer Connections.

Table 2.10. System Computer Connections

Part	Description
USB 2.0 Port	Connects a USB 2.0 cable to the USB 2.0 connector on the instrument's Left-Side Connectors
USB 3.0 Port	Connects a USB 3.0 HUB to the USB 3.0 HUB connector on the instrument's Left-Side Connectors

2.9 SPECIFICATIONS

This section lists the technical specifications for the Phenolmager HT instrument. Technical specifications are subject to change without notice.

NOTE *Phenolmager HT is for research use only. Not for use in diagnostic procedures.*

SYSTEM	
Size (H x W x L)	28" (72 cm) x 30" (77 cm) x 27.25" (69.2 5 cm)
Weight	185 lbs. (84 kg)
Spectral Range	440 nm to 780 nm
Objectives	4x, 10x, and 20x
Pixel Resolution	0.25 µm, 0.5 µm, or 1.0 µm
File Format	Akoya Biosciences. qptiff format for whole slide scans, proprietary .im3 file format for multispectral data; 24- bit Windows-compatible bitmap for RGB/Mono imagery
Operating System	Microsoft Windows 10, 64-bit
RAM	32 GB

ENVIRONMENTAL	
The location where the PhenolMager HT will be installed must meet the following specifications:	
Operating Temperature*	64°F - 83°F (18°C - 28°C)
Storage Temperature*	59°F - 86°F (15°C - 30°C)
Shipping Temperature** (up to 72 hours max)	14°F - 113°F (-10°C - +45°C)
Operating Humidity**	10% - 60% non-condensing
Storage Humidity**	10% - 60% relative humidity, non-condensing
Altitude	Up to 2000 meters (6561 feet)
Vibration	Up to 0.4mm/sec @ 8-80Hz
Pollution Degree	2

*Ambient temperature should not vary >2°C (4°F) over 24 hours.

**Do not expose the instrument to humidity >80% or temperature outside of -10°C to 45°C (14°F to 113°F) for more than 72 hours.

NOTE For indoor use only.

- **A clean, dust-free room is recommended.**
- **Do not operate the system in an environment with explosive or flammable gases.**
- **Do not operate the system in places where it may be splashed with liquid.**

ELECTRICAL	
One properly grounded AC power outlet for the computer, monitor, and instrument must be located within 6 feet (1.8 m) of the location.	
Input Voltage	100 to 230 VAC ($\pm 10\%$), 500 W, 50/60 Hz System does not have transient overvoltage protection.
Computer Interface	USB 2.0, USB 3.0
Fuse	4A Littelfuse [®] 250V, 5 mm x 20 mm

BARCODE READER LIGHT SOURCE (white laser diode)	
Wavelength	400 nm to 750 nm
Power Output	70 mAmps

Chapter 3: Operation

This chapter describes how to operate the PhenolMager HT hardware. It is important to read and understand the [Instrument Safety](#) section before using the system. The [Hardware Reference](#) chapter provides a description of each hardware component in the PhenolMager HT system.

Section	PAGE
3.1. System Startup	35
3.2. System Shutdown	36
3.3. Inspecting Slides and Slide Carriers	36
3.4. Loading Slides into the Slide Carriers	36
3.5. Loading Slide Carriers into the Slide Carrier Hotel	37
3.6. Removing Slide Carriers from the Slide Carrier Hotel	38
3.7. Removing Slides from the Slide Carriers	39


3.1 SYSTEM STARTUP

This section describes the procedures needed to start the PhenolMager HT system.

3.1.1 Turn on the PhenolMager HT Instrument

To start the PhenolMager HT instrument:

1. If necessary, plug the PhenolMager HT power cord into an appropriate power outlet.

	<p><i>Warning:</i></p> <ul style="list-style-type: none"> • Use only the power supply cord set provided with the PhenolMager HT system. If the correct cord set for the location was not provided, contact Akoya Biosciences Technical Support at support@akoyabio.com for a replacement. Do not use power supply cords with inadequate ratings. • Use only a properly grounded power outlet when connecting the system to power. • Appliance inlet is a disconnecting device. Place device or equipment in a manner so that disconnecting device is accessible at all times.
---	--

2. Turn on the computer and allow Windows 10 to start.
3. Switch the System Power Switch to the ON (I) position. The Status Progress Display flashes green to indicate a quiescent state. The Door Sensor is responsive and functional.

3.1.2 Launch the PhenolMager HT Software

To launch the PhenolMager HT software:

1. Double-click the PhenolMager HT icon on the Windows 10 desktop.
2. The PhenolMager HT home page opens. (See [4.1. Software Overview](#).)

3.2 SYSTEM SHUTDOWN

To shut down the Phenolmager HT system:

1. Exit the Phenolmager HT software. If open, the System Door closes.
2. Select Shut Down from the Windows Start Menu to power down the computer.
3. Switch the System Power Switch on the right side of the instrument to the OFF (O) position. The status lights turn off and the system shuts down.

NOTE For best performance, power cycle the Phenolmager HT system periodically (a few times per week) or if it has been idle for an extended period of time.

3.3 INSPECTING SLIDES AND SLIDE CARRIERS

Before inserting slides into slide carriers, both should be inspected for potential defects. This section describes the steps needed to properly inspect the slides and slide carriers.

3.3.1 Inspecting Slides

- Verify the slides meet the required formats and dimensions. (See [Table 2.9](#).)
- Do not use broken or damaged slides or slides with broken or damaged coverslips.
- Use only slides that are free of debris, fingerprints, and dust.

3.3.2 Inspecting Slide Carriers

- Verify the slide Carrier Tab Cover is secure.
- Do not use slide carriers that are warped or bent.
- If any sticky residue is on the carrier handle or outer surfaces, then clean the handle or surfaces before use. (See [5.4. Cleaning the Slide Carriers](#).)

3.4 LOADING SLIDES INTO THE SLIDE CARRIERS

After inspecting the Slides and Slide Carriers, load the slides into the Slide Carriers.

To load a slide into a Slide Carrier:

1. Place the Slide Carrier onto a flat surface.
2. Hold the microscope slide by the label end with the coverslip side facing up.
3. Gently push the opposite end of the slide into the preferred Slide Position until the slide is up against the Spring-Loaded Tab. (See [Figure 3.1](#).)




Figure 3.1. Loading a slide into a Slide Carrier.

4. Gently lay the label end of the slide completely into the slide position so the slide lays flat in the carrier. As the slide is released, the tab should push the slide against the top wall of the slide position.
5. After the slide is secured by the Spring-Loaded Tab, gently push the slide to the right so its long edge is flush with the right wall of the slide position.

NOTE *Slide indexing is critical to the PhenolMager HT slide detection process. Best practice recommendation is to index slides in the Slide Carrier to the right to ensure proper detection and scanning of all loaded slides.*

6. Load up to four slides into each Slide Carrier. Slides need not be contiguous.
7. Click the Enter Slide IDs button. See Slide IDs in the PhenolMager HT software. The Enter Slide IDs dialog opens.
8. Type the applicable slide IDs into the text boxes and click the OK button.
9. To load additional slides into another Slide Carrier, repeat steps 1-8 above. Slide Carriers may be stacked on top of one another for easy handling and storage. (See Figure 3.2.)



Figure 3.2. Slide Carrier stack.

3.5 LOADING SLIDE CARRIERS INTO THE SLIDE CARRIER HOTEL

The Slide Carrier Hotel stores Slide Carriers before and after microscope slide scanning. The Slide Carrier Hotel can store up to 20 Slide Carriers for a total of 80 slides. The Hotel Slots are numbered 1-20, starting at the top of the hotel.

NOTE *Slide Carriers can be loaded into the Slide Carrier Hotel before or after launching the software.*

To load a Slide Carrier into the Slide Carrier Hotel:

1. Inspect the Slide Carrier. (See [3.3.2. Inspecting Slide Carriers](#)). Notice the ends of the Slide Carrier. (See Figure 3.3.)
2. Verify the PhenolMager HT System is on. (See [3.1. System Startup](#))
3. If the System Door is closed, place a hand in front of the Door Sensor to open it. The System Door will slide to the left to reveal the Slide Carrier Hotel.
4. Hold the Slide Carrier by the Carrier Handle with the slide labels facing upwards into a Hotel Slot.


5. Identify which end of the Slide Carrier is marked with the Insert Indicator. ()
6. Gently insert the end of the Slide Carrier with the Insert Indicator into a Hotel slot.



Figure 3.3. Ends of a Slide Carrier.

7. When the Slide Carrier is halfway into the Hotel Slot, the carrier engages a roller and creates a small amount of friction. Push the carrier further into the slot until it is fully seated. (See Figure 3.4.) Hold the Slide Carrier by the Carrier Handle with the slide labels facing upwards into a Hotel Slot.




Figure 3.4. Pushing the Slide Carrier in until fully seated.

- If the carrier is inserted correctly, a click noise indicates that the carrier has been seated correctly into the slot. The Slide Carrier Status Light next to the associated Hotel Slot turns white.
 - If the carrier is not inserted correctly, the Slide Carrier Status Light next to the Hotel Slot blinks red until the Slide Carrier is removed from the slot or is positioned correctly.
8. Continue to load up to 20 Slide Carriers into the slots.
 9. When completed, place a hand in front of the Door Sensor to close the System Door.

3.6 REMOVING SLIDE CARRIERS FROM THE SLIDE CARRIER HOTEL

When slide scanning is complete, the Slide Carriers can be removed from the Slide Carrier Hotel.

To remove a Slide Carrier from the hotel:

1. Pull the Carrier Handle until the Slide Carrier is free from its Hotel Slot. The Slide Carrier Status Lights next to the Hotel Slot turns off, indicating no Slide Carrier is inside the Hotel Slot.
2. Lay the Slide Carrier onto a flat surface.

3.7 REMOVING SLIDES FROM THE SLIDE CARRIERS

After Slide Carriers are removed from the Slide Carrier Hotel, the microscopic slides can be removed from the Slide Carriers.

To remove a slide from a Slide Carrier:

1. Gently push against the slide label end of the microscope slide using an index finger, compressing the tab on the far end of the slide.
2. Using the same finger, lift the label end of the microscope slide from the slot.
3. Grab the label end of the slide with thumb and index finger.
4. Remove the remaining end of the slide from the Slide Carrier.

Chapter 4: Software Operation

This chapter describes how to operate the Phenolmager HT software.

Section	Page
4.1. Software Overview	40
4.2. System Dashboard	43
4.3. Creating and Editing Protocols	48
4.4. Scanning Slides	69



Figure 4.1. Phenolmager HT software.

4.1 SOFTWARE OVERVIEW

The Home Page links to the pages needed to maintain and run the Phenolmager HT. The Home Page is seen in Figure 4.2.




Figure 4.2. Phenolmager HT software home page.

The Home page buttons link to the following pages:

- **Check Dashboard** – Check the remaining disk space and acquire references.
- **Edit Protocol** – Create or edit brightfield and fluorescent protocols and studies.
- **Scan Slides** – Select the rules to scan the slides and scan the slides based on those rules.
- **Launch Phenochart** – Launch the Phenochart viewing software.

The Home pages can be accessed by using the Gear Menu located in the upper-right corner of the window. (See Figures 4.2 and 4.3.) See Section 4.1.1: Gear Menu for more information.




Figure 4.3. Gear menu.

Additionally, the Gear menu can be used to check the version number or perform other functions and launch the online help.

Below is an overview of each page.

Check Dashboard

Use the Dashboard to check the status of the system prior to imaging. It includes the following:

- Disk space—checks the available space on the disk drive where images will be stored.
- Acquire References—views and acquires brightfield references and fluorescent references.

For detailed information on how to use the System Dashboard, see [Section 4.2. System Dashboard](#).

Edit Protocol

Use the Edit Protocol page to define protocols and take snapshots.

A protocol describes how a sample is to be imaged, including the imaging mode (brightfield or fluorescence), and the spatial resolution (magnification) for the whole slide scan and for multispectral regions of interest. For fluorescent imaging, it also describes the exposure times and what bands to use for focusing and imaging.

For detailed information on how to create and edit protocols and studies, see [Section 4.3 Creating and Editing Protocols](#).

Scan Slides

Use the Scan Slides page to declare how slides should be imaged. Slide scanning can be started and stopped from this page.

For detailed information on how to perform whole slide scans, see [Section 4.4 Scanning Slides](#).

Launch Phenochart

Use the Launch Phenochart to launch the Phenochart program.

For detailed information on how to use Phenochart, see the [Phenochart User Guide](#) and online help from within the Phenochart program.

4.1.1 Gear menu

Use the Gear menu in the upper right-hand corner of the Home page to:

- Link to the following page:
 - System Dashboard
 - Edit Protocol
 - Scan Slides
 - Launch Phenochart
- Find out what's new in the current installed version.
- Launch the online Help system.
- Open the Preferences dialog.

- Open the Phenolmager HT software log for reference or troubleshooting.
- Open the About window to view the current software version.
- View contact information for Akoya Biosciences technical support.

4.2 SYSTEM DASHBOARD

Use the System Dashboard shown in Figure 4.4 to:

- See the available [Disk Space](#).
- View and acquire [Brightfield References](#).
- View and acquire [Fluorescence References](#).




Figure 4.4. System Dashboard window.

4.2.1 Disk Space

The disk space bar shows where the data is stored and how much space is available on that drive.

4.2.2 Brightfield References

Viewing Brightfield References

Click View to see reference imagery and information for the label image, brightfield overview, color image and multispectral image. (See Figure 4.5.)

- Export For Diagnostics will save this image to aid in technical support.
- Show raw images will display the images without scaling, which is useful for technical support.

- Compensation Information describes how these reference images are applied. (See [Section 4.2.4. Compensation Information](#))




Figure 4.5. View Brightfield References window.

Acquiring Brightfield References

Click Acquire to take new brightfield references. (See Figure 4.6.) Use the stage control to:

- Load the Reference Carrier onto the slide.
- Move around to a clean area in the live view (no coverslip lines, tissue, or label) using the stage navigation tool.
- Click within the slide to change positions. Refine the position by using arrow keys for small movements, and control + arrow for slightly larger movements.
- Click OK to take references.




Figure 4.6. Prepare for Brightfield Reference Acquisition.

4.2.3 Fluorescence References

Viewing Fluorescence References

Click View to see reference imagery and information for the overview, 20X and 40X resolution references for each filter. (See Figure 4.7.)

- Export For Diagnostics will save this image to aid in technical support.
- Compensation Information describes how these reference images are applied and the Fitted Flatness assessment. (See [4.2.4. Compensation Information](#))
- Responsivity Information describes how the system response is measured.




Figure 4.7. View Fluorescence References window.

Responsivity Information

Responsivity values recorded by Phenolmager allow the monitoring of the overall system response to a standard fluorescent artifact over time. Responsivity is affected by system components and alignments related to fluorescence illumination and detection, so a substantial decrease in Responsivity may indicate that the reference images need a closer look by the Akoya support team.

The measurements are made using acrylic fluorescent artifacts, which can vary in intensity between systems but do not vary significantly over time. Responsivity values are assessed for each combination of objective and filter band that can be used for fluorescence imaging each time Fluorescence References are acquired. For each band, several images are acquired of the fluorescent artifact, then a signal intensity is calculated by averaging the mean intensity within the central 1/16th area patch of the images. That signal intensity is then normalized by the exposure time and other camera settings to produce a Responsivity value with units of counts per second.

Acquiring Fluorescence References

Click Acquire to generate new fluorescent references. We recommend acquiring references for all filters at once but can take references for a single filter if necessary. Acquiring references for all filters can take over 45 minutes, and the system can be left unattended during this time.

Use the stage control to:

- Load the Reference Carrier onto the stage.
- Move to a tissue-filled location on the fluorescent reference slide using the stage navigation tool.

- Click within the slide to change positions. The position can be refined by using arrow keys for small movements, and control + arrow for slightly larger movements.
- Click OK to acquire references.




Figure 4.8. Acquire Fluorescence References window.

4.2.4 Compensation Information

Brightfield

Acquired images of a sample are normalized (divided) by the reference images on a pixel-by-pixel, wavelength-by-wavelength basis. This yields the sample transmission T , which is in the range 0–1. This transmission (T) is then mapped:

- For multispectral images, the Optical Density is calculated—this is based on $\text{Log}_{10}(T)$.
- For simple color images such as a whole slide scan, the transmission is mapped from 0-255.

Fluorescence Protocols

Acquired images are normalized by a shading pattern derived from the reference images to produce flat-fielded imagery. The shading pattern, which is the reference image divided by the mean intensity in the center, is applied on a pixel-by-pixel basis, with one pattern per epi-filter. Here, center means the central 1/16th of the image area. While the exact shape of the shading pattern varies per instrument, the overall effect is to increase the signal near the image edges, and to do little or nothing to the signal from the center of the image.

The PhenolMager provides a Fitted Flatness metric to estimate the quality of flat fielding after dividing by the reference image (metric shown in Figure 4.7). The Fitted Flatness metric is measured from overlapping images of the fluorescence reference slide acquired during fluorescence reference acquisition ([Section 4.2.3](#)), after applying shading compensation. The overlapping, compensated images are divided into small strips to compare the same area of tissue imaged in the center of the field vs the edges of the field. The Fitted Flatness metric reports the worst measured intensity difference observed when comparing the center of the field to any of the four edges of the field.

4.3 CREATING AND EDITING PROTOCOLS

4.3.1 Studies

Phenolmager HT stores scan data by Study.

A Study is a group of slides that belong together. This could be an experimental study (e.g., Ki67 markers in breast cancer tissue), all slides from one source, or other groupings. Each study contains one or more slides. (See Figure 4.9). Each slide may be scanned more than once, if needed.

- The default location for a study is D:\Data\PhenolmagerHT\[Study] (where [Study] is the name of the study).
- Whole-slide scans and supporting imagery acquired from specific slides are saved to slide-specific subfolders in the main study folder. See the section on Scanning Slides for more details about imagery.
- Multispectral Fields acquired from a particular slide scan will be stored in an MSI folder within the scan folder.

A protocol defines the set of rules to be used during whole slide and multispectral region acquisition, including imaging mode, pixel resolution, filter cubes, exposure times, and other parameters.

- Protocols have the file extension “.ppr” and are saved in D:\Data\PhenolmagerHT\[Study]. Figure 4.10 shows the QPTIFF format for the slides.




Figure 4.9. Study structure.



Figure 4.10. Slide QPTIFF.

4.3.2 Creating Protocols

Before creating a protocol, it is helpful to know how the slide was stained.

- Brightfield protocols are used to acquire imagery from slides stained with H&E or conventional chromogenic IHC methods.
- Fluorescence protocols are used with Akoya Biosciences PhenoCode Signature panels, and its Opal dyes.

To create a new protocol:

1. Select Edit Protocol from the PhenolMager HT Home Page.
2. Click the New button. The Create New Protocol window opens. (See Figure 4.11.)




Figure 4.11. Create New Protocol window.

3. Enter a Protocol name.
4. Under Imaging Mode, select either Brightfield or Fluorescence.
5. If the system has been upgraded for multispectral whole slide scans, choose the type of staining kit used for this protocol, or choose to select custom filters for imaging.
6. Select a previously created Study or create a new Study.

- To select a previously created study, click on the study in the Available Studies list. This will highlight the study.
 - To create a new study, enter the Study Name in the text box and click Create Study. The new study will be added to the Available Studies list where the study can be selected.
7. Click Create Protocol or Select Scan Bands to create the protocol in the selected study.
 8. If choosing to select custom filters, the Select Scan Bands window opens. Select the bands to use for imaging and click OK.

To load an existing protocol, click Load and select the protocol from the study folder.

4.3.3 Editing Protocols

After the brightfield or fluorescence protocol has been created and assigned to a study, use the Edit Protocol screen to add specific details to the protocol. The next sections are organized by two different types of protocols: brightfield and fluorescence.

Brightfield Protocols

After the brightfield protocol has been created and assigned to a study, the Edit Protocol window (for brightfield protocols) opens. (See Figure 4.12.)




Figure 4.12. Edit Protocol Window.

Under Slide Scan Settings, select Pixel Resolution in the drop-down list to select the resolution to apply to the whole slide. (See Figure 4.13.) Choose 1 µm per pixel (nominally 10x), 0.5 µm per pixel (nominally 20x), or 0.25 µm per pixel (nominally 40x).




Figure 4.13. Slide Scan Settings - Pixel Resolution drop-down list.

Under Multispectral Fields, select the Pixel Resolution in the drop-down list to use for imaging. (See Figure 4.14.) Again, available resolutions are 1 µm per pixel (nominally 10x), 0.5 µm per pixel (nominally 20x), or 0.25 µm per pixel (nominally 40x).




Figure 4.14. Multispectral Field Settings - Pixel Resolution drop-down list.

Advanced Settings

There may be situations when fine-tuning how PhenolMager HT scans the slides is necessary. Figure 4.15 showcases the Advanced Settings button. Click the Advanced Settings button to adjust the following functions:

- **Coverslip Thickness** - PhenolMager HT is designed for use with standard #1.5 coverslips, where the distance between the tissue and the top of the sample is approximately 170 microns. While this is the optimal setup, PhenolMager HT can handle other scenarios.
 - If using thinner coverslips (#0), select 100 µm (#0 or similar) for a Coverslip Thickness.
 - If using thicker coverslips or using a significant amount of mounting media between the sample and coverslip, select the option that best matches the slides.

NOTE *PhenolMager HT is capable of handling slides with up to 420 microns of combined coverslip or mounting medium above the sample, though imagery at these extremes may be suboptimal.*

- **Scan within the entire coverslip region** - If PhenolMager HT cannot find the entire tissue sample, check the Scan within the entire coverslip region checkbox. This will increase scan time and file size but will enable complete scanning of difficult samples.
- **Sample is a TMA** - PhenolMager HT has a specialized algorithm for finding TMA samples. If samples are TMAs, select Sample is a TMA.
- **Image Compression** - By default, brightfield whole slide scans in PhenolMager HT are .jpg compressed to save disk space. Adjust the image quality as needed; higher quality will result in larger files. If less compression is preferred, select LZW rather than JPEG.
- **Label Barcode Reading** - PhenolMager HT can decode the barcode on the slide label and save it as part of the scan file. To do this, select the type of barcode being used. By default, PhenolMager HT will not attempt to decode a barcode unless this option is specifically selected.
 - If only 1D barcodes, select All 1D barcodes (autodetect type).
 - If there are a variety of 2D barcodes in use, select All 2D barcodes (autodetect type).

NOTE *If the lab only uses a specific type of 2D barcodes, the decoding time can be decreased by selecting the specific type. PhenolMager HT supports decoding of Data Matrix, PDF 417, QR Code, Aztec, and Maxicode.*

- If using a combination of 1D and 2D barcodes, select All barcodes (autodetect type). This option will increase scan time and file size and should only be used if necessary.




Figure 4.15. Advanced Settings window.

Click Save to save the protocol. (See Figure 4.16.) The study previously chosen will be automatically selected. The protocol name or study can be changed if needed.




Figure 4.16. Save Protocol window

Brightfield Snapshots

Once the scan and multispectral field resolutions have been selected, example snapshots of the slides can be taken.

Click Take Snapshots to load the Brightfield Snapshots editor. (See Figure 4.17.)




Figure 4.17. Brightfield Snapshots.

If a carrier is not already on stage, click Load Carrier to select a carrier that contains the slides to be imaged. (See Figure 4.18.)




Figure 4.18. Select Carrier.

Use the navigation tool to select an area on the slide. (See Figure 4.19.) Click on the slide to change positions.




Figure 4.19. Navigation Tool.

For easier navigation, click Take Overview to acquire imagery of all slides in the carrier. An option to do this automatically can be selected in the Preferences dialog in the Gear menu. PhenolMager HT will then automatically acquire overview imagery whenever a carrier is loaded.

Click on the tissue to change positions.




Figure 4.20. Navigation tool with Overview.

Refine the position by using arrow keys (for small movements) or control + arrow (for slightly larger movements).

Click Autofocus or use the Stage Height slider to bring the live view into focus.

Click on a feature within the live view to center on that location or click and drag within the live view to change locations.

To take a snapshot of the current live view, select either the whole slide or the MSI regions radio button (see Figure 4.21), pick a base file name, and click Snap.




Figure 4.21. Whole Slide and MSI Regions radio buttons.

Navigate to new places, focus as needed, and take as many snapshots as required. Snapshots will be numbered incrementally.

Click Back when ready to return to the Edit Protocol window.

Fluorescence Protocols

When creating a fluorescence protocol, slides are required that include positive expression in all markers of interest to set suitable exposure times.

After creating the fluorescence protocol and assigning it to a study, the Edit Protocol window (for fluorescence protocols) opens. (See Figure 4.22.)




Figure 4.22. Fluorescence Protocol editor.

Under Overview Scan Rules, select the desired filter that will be used to help find tissue on the slide. This will typically be DAPI counterstain. If there is no counterstain, then choose a filter that aligns with the expression of the most common stain or auto-fluorescence.



Figure 4.23. Overview Scan Rules.

Under Slide Scan Settings, choose the Pixel Resolution to be used while imaging the slide. (Figure 4.23). Options include:

- 1.0 µm per pixel (nominally 10X).
- 0.5 µm per pixel (nominally 20X).
- 0.5 µm per pixel (nominally 20X, binned from a 40X acquisition).
- 0.25 µm per pixel (nominally 40X).

Choosing higher resolution scans will result in larger file sizes and longer scan times.




Figure 4.24. Slide Scan settings.

NOTE Each time the pixel resolution is changed the exposure times will need to be reset, as displayed in Figure 4.24.

Exposure Times




Figure 4.25. Exposure times require update window.

Functions on the Exposure Times window are described below. (See Table 4.1.)

Table 4.1. Exposure Times window functionality

Button	Description
Auto-Update	Option to auto-update the exposures. This is recommended as a first step to give a close approximation of the exposure times. We then recommend using the exposures editor to fine-tune the exposures.
Manual Update	Option to manually update the exposures. In this case, the current exposures remain unchanged and must use the exposures editor to obtain valid exposures.
Cancel	Option to cancel and leave the pixel resolution unchanged.

If planning to take Multispectral Fields, choose the Pixel Resolution that should be used for fields imaging. Available resolutions are 10x, 20x, 20x binned, or 40x.




Figure 4.26. Multispectral Field Settings.

The Select Scan Bands button provides the option to choose which bands used to take the slide scan. Imagery for each band will be taken in the order it is shown. By default, PhenolMager HT will take a slide scan using the best selection available for the staining kit, or the bands selected for a custom scan.

To remove bands where fluorescent signals are not expressing, uncheck the band. Drag-and-drop to change the order in which they are acquired. The colors associated with each band can be changed; these are the colors that Phenochart will use to display the scan.

NOTE *If the Overview or Focus bands are removed, PhenolMager HT will select new bands. Visit Advanced Settings to confirm these selections.*




Figure 4.27. Select scan bands.

Limit the filters displayed in the Select Scan Bands dialog (Figure 4.27).

- Show Akoya Scan Filters to limit the available filters to installed Akoya filters that are appropriate for slide scanning.
- Show All Akoya Filters to display all installed Akoya filters, both for scanning and taking multispectral fields (Figure 4.28.)
- Show All Filters to display all installed filters, including custom filters that have been installed.




Figure 4.28. Select scan bands—show drop-down list.

The Select Field Bands button provides the option to choose which bands are used to take the multispectral fields. Akoya Biosciences strongly recommends maintaining the default band configuration. Non-signal bands are integral to multispectral stain extraction and unmixing.

To add or remove bands, select or de-select the filter. If there is more than one band available for a filter, select the band from the drop-down list. Drag-and-drop to change the order in which they are acquired. (See Figurer 4.29.)

NOTE *If the Focus band is removed, Phenolmager HT will select a new band. Visit Advanced Settings to confirm this selection.*




Figure 4.29. Select Field Bands – band selection.

Limit the filters displayed in the Select Field Bands dialog.

- Show Akoya MSI Filters to limit the available filters to installed Akoya filters that are appropriate for multispectral fields. (See Figure 4.30.)
- Show All Akoya Filters to display all installed Akoya filters, both for scanning and taking multispectral fields.
- Show All Filters to display all installed filters, including custom filters that have been installed.




Figure 4.30. Select Field Bands – MSI filters.

Data Format and Unmixing

PhenolMager HT can save fluorescent scan data in several formats. By default, it uses “extended range” which provides a safety margin against accidental over-exposure. (See [1.3.6 Saturation and Dynamic Range](#).) The Standard format from Version 1 software may be used to save disk space. Finally, PhenolMager HT can spectrally unmix the data as part of the scan operation and record the marker names in the dataset.

To select between these modes, click the Data Format / Unmixing button on the Fluorescence Protocol Editor. This will bring up a dialog to select the data format, and an option to specify unmixing. (See Figure 4.31.)




Figure 4.31. Extended Range selection.

This shows the default selection and posts a reminder that this format uses more disk space than Standard. To select between the data formats, use the drop-down tool. (See Figure 4.32.)




Figure 4.32. Select Data Format and Unmixing.

When Standard is selected, there is no reminder about disk space usage. (See Figure 4.33.)




Figure 4.33. Standard Format Selection.

When Unmixed Data format is selected, a similar warning is shown, along with a notice to load an unmixing library. (Figure 4.34.) The unmixing library is contained in a .umxlib file that must be created in inForm and then exported for use by PhenolMager HT. It describes the properties of the dyes and sample autofluorescence that will be used to unmix the sample.




Figure 4.34. Unmixing Selection.

Click the Load Unmixing Library button and select a .umxlib file to use for the protocol. Once it has been loaded, there will be a prompt to fill in the names of the markers associated with each dye channel.

NOTE This enables the same .umxlib to be used for a variety of protocols, targeting different markers. (See Figure 4.35.)




Figure 4.35. Band name assignment.

While the example above shows a protocol for imaging 6 markers plus DAPI and sample AF, the protocol may contain fewer markers if it has DAPI and sample AF and at least 1 marker channel; and the unmixing library has entries for all the channels used.

This enables the creation of an unmixing library once, with all channels, that can be used on a range of protocols with fewer channels when imaging samples with fewer markers.

The Fluorescent Protocol Editor shows the data format, with information about the unmixing library and marker names when unmixed data format is selected. (See Figure 4.36.)




Figure 4.36. Fluorescent Protocol Editor.

Advanced Settings

There may be situations when fine-tuning how PhenolMager HT scans the tissue is necessary. Click the Advanced Settings button to adjust the following functions:

- PhenolMager HT is designed for use with standard #1.5 coverslips, where the distance between the tissue and the top of the sample is approximately 170 microns. While this is the optimal setup, PhenolMager HT can handle other scenarios.
 - If using thinner coverslips (#0), select “100 um (#0 or similar)” for a Coverslip Thickness.
 - If using thicker coverslips or have a significant amount of mounting media between the sample and coverslip, select the option that best matches the slides.

NOTE *PhenolMager HT is capable of handling slides with up to 420 microns of combined coverslip or mounting medium above the sample, though imagery at these extremes may be suboptimal.*

- Specify which filter will be used for setting focus during whole slide scans and for multispectral field imaging. Choose the filters where the counterstain expresses. This is normally DAPI.
- Saturation protection prevents overexposure when acquiring multispectral field imagery. Akoya Biosciences highly discourages deactivating the Use Saturation Protection setting.

- If Phenolmager HT cannot find the entire tissue sample, check the Scan within the entire coverslip region checkbox. This will increase scan time and file size but will enable complete scanning of difficult samples.

NOTE *If there is Sharpie markup on the fluorescent slide, Phenolmager HT will use that as the region instead of the coverslip.*

- Phenolmager HT has a specialized algorithm for finding TMA samples. If the samples are TMAs, select Sample is a TMA.
- Label Barcode Reading
 - Phenolmager HT can decode the barcode on the slide label and save it as part of the scan file, select the type of barcode being used. By default, Phenolmager HT will not attempt to decode a barcode unless this option is specifically selected.
 - If only 1D barcodes are used, select All 1D barcodes (autodetect type.)
 - If a variety of 2D barcodes is being used, select All 2D barcodes (autodetect type.)

NOTE *If the lab only uses a specific type of 2D barcodes, decrease the decoding time by selecting that specific type. Phenolmager HT supports decoding of Data Matrix, PDF 417, QR Code, Aztec, and Maxicode.*

- If a combination of 1D and 2D barcodes is used, select All barcodes (autodetect type). This option will increase scan time and file size and should only be used if necessary. (See Figure 4.37.)




Figure 4.37. Fluorescence—Advanced Settings.

After selecting the settings, use the Exposures Editor to set exposures for the protocol.

4.3.4 Setting Exposures

After Selecting the scan and multispectral field resolutions, filters, and bands, exposures need to be set. Click Set Scan Exposures or Set Field Exposures in the Fluorescence Protocol Editor to set exposure times for scanning and field acquisition. Figure 4.38 shows the Edit Exposures window.




Figure 4.38. Edit Exposures.

If a carrier is not already on stage, click Load Carrier to select a carrier that contains representative slides. Select the carrier from the list and load it to the stage. (See Figure 4.39.)




Figure 4.39. Select Carrier.

Use the navigation tool to select a slide. A live view of the first filter for the imaging rules will be shown. In most cases, this will be DAPI. The highlighted entry in the table identifies which filter or band is currently shown. (Figure 4.40.)




Figure 4.40. Navigation tool.

For easier navigation, click Take Overview to acquire digitally enhanced imagery of all slides in the carrier. (See Figure 4.41.) That option in the Preferences dialog can be selected in the Gear menu. Phenolmager HT will then automatically acquire overview imagery whenever a carrier is loaded.




Figure 4.41. Navigation tool with overview.

Click on the tissue to change positions.

Refine the position by using arrow keys (for small movements) or Ctrl + arrow (for slightly larger movements).

Click on a feature within the live view to center on that location or click-and-drag within the live view to change locations.

Click Autofocus or use the Stage Height slider to bring the live view into focus. Click Autoexpose to have the system find the best exposure for that filter/band.

- After auto-exposing, auto-focus and auto-expose again to refine the focus and exposure estimates.
- Override the auto-exposure value by typing a value in the highlighted cell. Values must be between 0.1 and 2000 ms.
- To turn off the fluorescence illumination and live view, click Close Shutter. Then re-open the shutter to see a live view and set exposures.

If this protocol is intended for multispectral field acquisition, set the multispectral field exposures at this time. In most cases, the same area of the slide can be used to set the multispectral field exposures. Select the Set Field Exposures radio button to set exposures for

multispectral fields. Focus and click Autoexpose.

- As with whole slide scan exposures, after auto-exposing, auto-focus and autoexpose again to refine the focus and exposure estimates.
- Override the auto-exposure value by typing a value in the highlighted cell. Values must be between 0.1 and 2000 ms.
- To turn off the fluorescence illumination and live view, click Close Shutter. Re-open the shutter to see a live view and set exposures.

Repeat the steps above for all filters and bands in the protocol. Locations and/or slides may need to change to find the best signals for setting the exposures.

Filter Cube Recommendations when using Akoya Biosciences Opal Reagents

If using the PhenolMager HT system with Akoya Biosciences' Opal fluorophores (either as part of PhenoCode Signature Panels or as individual fluorophores), please refer to Table 4.2 and Table 4.3 for exposure time recommendations for acquiring spectral library images.

Table 4.2. Multispectral Slide Scan Band

Columns represent filter bands; rows represent Opal fluorophores.

		FILTER BANDS							
		DAPI	Opal 480	Opal 520	Opal 570	Opal 620	Opal 690	Opal 780	Sample AF
OPAL FLUOROPHORES	DAPI	DAPI							
	Opal 480		Opal 480						
	Opal 520			Opal 520					
	Opal 570				Opal 570				
	Opal 620					Opal 620			
	Opal 690						Opal 690		
	Opal 780							Opal 780	

When acquiring multispectral fields:

Table 4.3. Multispectral Field Band

Columns represent filter bands; rows represent Opal fluorophores.

OPAL FLUOROPHORES	FILTER BANDS								
		DAPI	Opal 480	FITC	Cy3	Texas Red	Cy5	Opal 780	
	DAPI	Blue	Dark Gray	Dark Gray	Light Gray	Light Gray	Light Gray	Light Gray	Light Gray
	Opal 480	Dark Gray	Cyan	Green	Light Gray	Light Gray	Light Gray	Light Gray	
	Opal 520	Dark Gray	Light Gray	Light Gray	Light Gray	Light Gray	Light Gray	Light Gray	
	Opal 540	Dark Gray	Cyan	Green	Yellow	Light Gray	Light Gray	Light Gray	
	Opal 570	Dark Gray	Dark Gray	Green	Yellow	Orange	Light Gray	Light Gray	
	Opal 620	Dark Gray	Dark Gray	Dark Gray	Yellow	Orange	Red	Light Gray	
	Opal 650	Dark Gray	Dark Gray	Dark Gray	Light Gray	Orange	Red	Light Gray	
	Opal 670	Dark Gray	Dark Gray	Dark Gray	Light Gray	Orange	Red	Light Gray	
Opal 690	Dark Gray	Dark Gray	Dark Gray	Light Gray	Light Gray	Red	Light Gray		
Opal 780	Dark Gray	Dark Gray	Dark Gray	Light Gray	Light Gray	Light Gray	Magenta		

Colors are signal bands: autoexpose or manually set an appropriate exposure per sample.

Dark Gray is an autofluorescence band: autoexpose or manually set an appropriate exposure per sample.

Light Gray is a non-signal band: manually set exposure according to Table 4.4.

Table 4.4. Exposure Times For Non-Signal Bands

	Resolution			
	10X	20X	20X (40X binned)	40X
Multispectral Slide Scans	6 ms	25 ms	2.5 ms	8 ms
Multispectral Fields	35 ms	150 ms	15 ms	50 ms

Snapshots

Once the exposures have been set, take snapshots to see sample imagery of how the exposures perform on a given area of the slide. Select either the Set Scan Exposures or the Set Field Exposures radio button, pick a base file name, and click the Snap button. (See Figure 4.42). Navigate to new places and take as many snapshots as desired. They will be numbered incrementally.




Figure 4.42. FL exposures panel.

Click the Back button when ready to return to the Edit Protocol window.

NOTE *The exposures just set are only valid for the resolutions that were selected. If the scan or multispectral field resolution is changed, then revisit the exposures editor to update the exposures.*

4.4 SCANNING SLIDES

4.4.1 The Carrier

Slides are loaded into PhenolMager HT via carriers, which are stacked into the Slide Carrier Hotel. (See Figure 4.43.)



Figure 4.43. Slide Carrier.

The hotel can hold up to 20 carriers, and each carrier can hold up to 4 slides. This means Phenolmager HT can be loaded with up to 80 slides at any given time. Those 80 slides can all be scanned based on the same set of rules, or the scan rules can be tailored to each slide and/or each carrier as needed.

4.4.2 Carrier Status

Within the software, each hotel slot is identified as Slot N, where N is 1-20 running top to bottom (Figure 4.44). When a Slide Carrier is loaded into the hotel, an LED light changes color to identify the status of that carrier.




Figure 4.44. Scan slides—carrier status.

The User Interface will reflect the current state of each slot in the hotel. If there is no carrier present, a striped display will result. (See Figure 4.45.)



Figure 4.45. Striped display indicating no carrier.

Carrier presence in the hotel is represented as indicated below in Figure 4.46.



Figure 4.46. Slide carrier indicators in the hotel.

Each carrier is represented by a carrier status icon, 4 slide status icons, and 4 slide IDs.

The carrier status icon corresponds to the hotel LED for that particular carrier. Each possible status is represented by a different color.

When the carrier is inserted into the hotel, a white icon indicates that the carrier is present but has no rules to use for slide scanning. If the icon is blinking red, the carrier has not been completely inserted; try pushing it in further.

Once the carrier has been inserted in the hotel, define the rules for the scan. The icon will change to blue when there is at least one slide ready to be scanned. Carriers that have been prioritized will be blinking blue.

If no slides on the carrier can be scanned with the current rules, the icon will change to yellow. This carrier will be ignored until the problem is resolved.








When scanning has started, the icon for the carrier on stage will be blinking green. When the carrier is returned to the hotel, the icon will either be green (for a successful scan) or orange (if there was a problem with the scan).

If the icon is solid red, that hotel slot's sensor is malfunctioning. Do not use this slot.

4.4.3 Slide Status

Each carrier also has a status icon for the 4 slides it contains. Table 4.5 shows the status icons and their corresponding descriptions for each slide.

Table 4.5. Slide status icons and their meanings

Status Icon	Description
	The slide has passed the requirements and is ready to be scanned.
	The slide was previously skipped due to insufficient disk space and is ready to be scanned. Alternately, there may be a warning related to instrument configuration. Hover over the icon to view details.
	<Blinking> The slide is currently being scanned.
	The slide was scanned successfully.
	The slide has rules, but they cannot be used. (They may be incomplete, in which case the rules will need to be edited for that slide.)
	The slide failed to scan correctly. Hover over the icon to view details of the error.
	The slide is being ignored. This is useful if the carrier has <4 slides or if that slide will not be processed at the current time.

4.4.4 Setting Up Scan Rules

Scan rules must be defined for each slide to process. Scan rules are complete when a study, protocol, task, and slide ID have been assigned.

Slide IDs

The ID must be entered manually for each slide in the carrier. The IDs may be entered before the carrier is placed into the hotel or enter IDs after the carrier has been placed into the hotel.

To enter the IDs before the carrier is placed into the hotel, click Enter Slide IDs. (See Figure 4.47.)




Figure 4.47. Enter Slide IDs button.

An Enter Slide ID window appears. (See Figure 4.48.)




Figure 4.48. Enter Slide IDs window.

Enter between 1-4 slide IDs for the slides and click OK. Immediately put the carrier in the hotel, and Phenolmager HT will assign the slide IDs that were entered to that carrier.

If the carrier is already present in the hotel, edit the carrier to add the slide IDs along with the other scanning rules.

Editing a Single Slide

To edit the scanning rules for a single slide, click the status icon for that slide. It will bring up the Edit Slides dialog box. (See Figure 4.49.)




Figure 4.49. Edit Slides.

For this slide, enter:

- **Task** - Choose Scan Slides or Acquire Fields.

NOTE Select Ignore to skip processing this particular slide.

- **Study** - Select a study.
- **Protocol** - Select any protocol saved within the study.
- **Slide ID** - Add the slide ID, if necessary.

To prioritize this slide's carrier to be scanned at the front of the line, check Prioritize the Containing Carrier.

Editing All Slides Within a Single Carrier

To edit the scanning rules for a single carrier, click the status icon for that carrier. The Edit Slides editor window is seen in Figure 4.50.




Figure 4.50. Edit Slides editor window.

For each slide, enter:

- **Task** - Choose Scan Slides, Acquire Fields, or Ignore. Select Ignore if there is no slide in the slot.
- **Study** - Select a study.
- **Protocol** - Select any protocol saved within the study.
- **Slide ID** - Add the slide ID, if necessary.

Check Prioritize This Carrier if the carrier should be scanned first.

Editing Multiple Carriers

To edit multiple carriers, click Configure Tasks. (See Figure 4.51.)




Figure 4.51. Configure Tasks button.

Select the carriers to be edited. Multiple selection is available. (See Figure 4.52.)




Figure 4.52. Configure Carriers.

Select how to process the slides.

To process all slides using the same rules, choose Using the Same Rules for all Slides and click Next. The Edit Slides window opens. (Figure 4.53.)




Figure 4.53. Edit Slides.

- Select the task, study, and protocol.
- To prioritize the selected carriers to be scanned at the front of the line, check Prioritize the Containing Carriers.

NOTE *The Slide ID cannot be edited because rules are applied across multiple slides. If slide IDs must be entered, edit the individual carriers.*

To use different processing rules for each carrier position, choose Using Different Rules for Each Carrier Slot and click Next For each position in the carrier, select the Task, Study, and Protocol. Prioritize if needed. Again, the Slide IDs cannot be edited because rules are applied across multiple slides. (See Figure 4.54.)




Figure 4.54. Edit Slides.

4.4.5 Assessing Disk Space Before Scanning

After setting up rules for all slides, assess whether any slides might be skipped due to insufficient disk space by viewing the disk space indicators at the bottom of the window. (See Figure 4.55.) Total available disk space is displayed at the bottom right. The bottom left text estimates how many slides will complete successfully given the available disk space.




Figure 4.55. Disk space indicators.

Since file sizes depend on the imaged area of each slide, which is not determined until a slide is loaded on the stage, predictions are shown for both typical scan areas (15 mm x 15 mm) and maximum scan areas (29 mm x 57 mm). If unsure of the scan areas, assume the maximum scan area to avoid skipped slides.

The text at the bottom of the screen will be red if any slides are predicted to be skipped even if all are typical size. The text will be blue if slides are only predicted to be skipped if all are maximum size. (See Figure 4.56.) Clear disk space to reduce the estimated number of skipped slides.

All 6 slides will complete if typical size with 202 GB extra.
 Data space warning: only 5 of 6 slides will complete if maximum size. Available disk space: 228 GB free.

Figure 4.56. Blue text indicating slides are only predicted to be skipped if all are maximum size.

4.4.6 Saving and Restoring Scan Setup

After setting up all the details for the slides, save work for later use. This may be useful if the scanning process is interrupted before the task is finished or if returning later to re-scan certain slides.

Save the slide IDs, studies, protocols, and tasks by clicking the Save Setup button. Phenolmager HT will also prompt to save before scanning (after clicking Scan). At this point, enter a name and some notes about this scan setup. (See Figure 4.57).




Figure 4.57. Scan Slides.

Saving a setup is not required, but it can help locate the information later if needed. (See Figure 4.58.)

- By default, this information will be deleted after 30 days, but can override the default and opt to keep the information until manually deleted.




Figure 4.58. Save Scan setup.

To restore an earlier scan setup, click Load Setup.

All available setups will be shown. Filter by date to show setups only from a specific time period. Click the headers to sort by date, name, notes, studies, etc.

NOTE *Entries can be deleted that are no longer needed by selecting those items and clicking Delete.*

To use a setup, select it from the dialog and click Load. (See Figure 4.59.)




Figure 4.59. Load Scan setup.

After selecting a setup, configure the carriers in PhenolMager HT to match the information in the setup.

If the carriers and slides are the same location as what was saved during setup, then click OK because everything is in the right place.

In other cases (for example, when putting the slides back in), this dialog will assist in figuring out where the slides belong.

1. When scan information is available and a carrier is present in the hotel, the system will show the slide IDs. (See Figure 4.60, Slot 1.)
2. When scan information is available but there is no carrier present in the hotel slot, the slot will be labeled as Missing in red. (See Figure 4.60, Slot 5.) The slide IDs for that carrier are listed, so one can locate the correct slides, place them in the carrier, and insert the carrier into the hotel.
3. When no information is available for a carrier that is present in the hotel, the system will show italicized Slide N text for the slide IDs (e.g., Slots 7 and 9). Slides in these carriers were not present in the original setup.

For every item in the scan setup (meaning anywhere the system has information on a carrier that should be present [Items 1 and 2]), before the carrier can be loaded, one of the following must be specified:

- **Use as saved:** The system will load all slide information including the scan task.
- **Ignore task:** The system will load all slide information except for the scan task. The scan task will be set to ignore. This is useful for loading the setup to rescan a few select items, or to change the task from how it was originally saved.
- **Don't restore:** No information for this carrier will be reloaded. This is useful if the carrier is not needed or is missing.




Figure 4.60. Slide setup.

The red text at the bottom of the screen will list any issues that need to be fixed before scanning.

Click OK when the selections have been finalized. The scan slides editor will be cleared and replaced with the selected scan information.

4.4.7 Scanning

The Scan button is enabled when at least one carrier is ready to be scanned. Click Scan to start scanning. If any carriers have been prioritized, these carriers will be processed first.

PhenolImager HT will report scanning progress for each slide on the progress dialog and on the front panel LEDs. (See Figure 4.61.)




Figure 4.61. Processing Slides window.

Click Pause Scanning to remove completed carriers and add new carriers to be scanned. When a scan is paused, Phenolmager HT will return the current carrier to the hotel. It will then be safe to open the door. Remove carriers that are completed and replace them with unscanned carriers. Edit those carriers to set them up for scanning and reprioritize carriers if needed. When ready, click Scan to resume scanning.

NOTE *If scanning using the Unmixed Data format (see Section 4.3.3), the last unmixed slide may still be processing after the final carrier is unloaded. In this case, the progress dialog will indicate status as “Performing unmixing on the final scan.” This step could take up to 120 minutes if using 40x scan resolution and large scan area.*

Chapter 5: Maintenance

This chapter includes procedures that are to be performed periodically, either to make the PhenolMager HT run better or to protect its components from damage.

Section	Page
5.1 Power Cycling the Instrument	81
5.2 Cleaning the Instrument Exterior	81
5.3. Cleaning the Monitor	81
5.4. Cleaning the Power and Communication Ports	82
5.5. Cleaning the Slide Carriers	82
5.6. Replacing the Fuses	82

5.1 POWER CYCLING THE INSTRUMENT



For best performance, power cycle the PhenolMager HT system periodically (a few times per week) or if it has been idle for an extended time. See section 3.1 and 3.2 for startup and shutdown procedures.

5.2 CLEANING THE INSTRUMENT EXTERIOR

Clean the PhenolMager HT instrument exterior as necessary.

The non-electrical exterior parts of the PhenolMager HT can be wiped down with a soft cloth using standard laboratory grade cleaning solutions including:



- 70% ethanol
- 10% bleach
- Clidox®
- Sporicidin®

Warning: DO NOT spray cleaning solutions directly onto the PhenolMager HT instrument. Sprays and liquids that come into contact with the PhenolMager HT instrument may result in damage to the system or electric shock. Always spray the cleaning solution onto a cloth and then wipe the instrument exterior with the cloth.

5.3 CLEANING THE MONITOR

Clean the monitor as necessary with a soft, lint-free cloth. If needed, dampen the cloth with water or an eyeglass cleaner.

Warning: DO NOT spray cleaning solutions directly onto the monitor screen. Sprays and liquids that come into contact with the screen may result in damage to the system or electric shock. Always spray the cleaning solution onto a cloth and then wipe the screens with the cloth.

5.4 CLEANING THE POWER AND COMMUNICATION PORTS

Dust can be removed from the power and communication ports with a gentle air stream from compressed air cans. DO NOT shake the can before use. Hold can at a distance from the ports to prevent condensation.



Warning: Turn off the electrical power to the Phenolmager HT system by shutting down the system (See [3.2. System Shutdown](#)) before cleaning any part of the instrument where electrical or fiber optic cables make connections.

5.5 CLEANING THE SLIDE CARRIERS

It is important to keep carriers free from sticky debris caused by slide labels, tape, or mounting media.

Slide carriers can be cleaned using soap and water with a soft cloth or using an ultrasonic bath.

5.6 REPLACING THE FUSES

For blown fuses and to order replacement fuses, contact Technical Support at support@akoyabio.com.



Warning: Electric Shock Hazard. DO NOT attempt to replace the fuses. Only qualified Akoya Biosciences service personnel can replace the fuses.

Appendix A: Akoya QPTIFF Specification

Background

This describes a TIFF format that Akoya Biosciences uses for its tissue images. The imagery may be a simple RGB image, a set of components extracted from a field or region, or a whole-slide scan. In the latter case, it may be a brightfield (BF) color RGB image or a multiband fluorescence (FL) image.

The goal is to use the same syntax and metadata for all these kinds of images, and minimize the semantic distinctions where possible. Specifically, an extracted component represents signal attributed to a stain or fluorescent dye in a sample. The signal values incorporate image preprocessing such as normalization for exposure time and spectral unmixing when that technique is used.

Data Format

The files are TIFF or BigTIFF images, depending on image size, with multiple images per file.

For images larger than about 2K x 2K pixels, tiled format is used, and the image is provided in several resolutions (pyramidal tiled images). Tile size is 512 x 512 pixels. Images smaller than 2K x 2K are stored in strips rather than in tiles.

The highest resolution (baseline) image(s) appear first in the file. For each resolution there are N baseline images where N depends on the contents. For BF images, N=1 and each image is an RGB image. For FL images or unmixed component images, N = number of bands, which is usually > 1, and each image is a grayscale image. The sample format for grayscale raw images can be unsigned 8-bit or unsigned 16-bit. For grayscale unmixed images, the sample format can be unsigned 16-bit, signed 16-bit, or 32-bit float.

A thumbnail RGB image is provided, and this is a good image to use as an icon in graphical image lists. Note that the thumbnail can be used for qualitative sample identification, but should not be used to assess relative brightness of different images. This comes after the baseline images, meaning it is the second image in BF (RGB) images, and the (N+1)st image for FL images or unmixed component images.

Next come the reduced-resolution images (if present). The pyramid contains enough levels that the image size is no larger than 2K x 2K at the coarsest resolution.

For whole-slide scans, there are two more non-tiled images after these: an optional RGB image of the label, and a macro (low-resolution) RGB image of the whole slide.

Overall, the arrangement is:

Table A.1. Images

Description	RGB/mono	Tile/Strip	Resolution	Notes
Baseline image	Varies	Varies	Full	Tiled if > 2K x 2K RGB for BF, else mono.
More full-resolution images	Mono	Varies	Full	If N > 1.
Thumbnail	RGB	Strip	~500 x 500	
Half-resolution images	Varies	Varies	Half	Only if baseline is tiled.
Quarter, eighth, etc.	Varies	Varies	Quarter, eighth, etc.	Continues until 2K x 2K or smaller
Macro (overview) image of whole slide	RGB	Strip	~2000 x 4000	Required for whole-slide scans. Optional for simple RGB images and extracted components.
Label image	RGB	Strip	~500 x 500	Optional, whole-slide scans.

Detection

Readers can recognize Akoya tissue images via the contents of the “Software” TIFF tag (see below). The file suffix is .qptiff for whole slide scans. Some TIFF tags contain “PerkinElmer” for historical reasons; they remain to retain 3rd party reader compatibility.

Metadata

Metadata is contained in two locations: within standard TIFF tags as listed in the table below, and within the ImageDescription string, using a set of XML tags described below. These are provided for each image (IFD) in the file, and describe that image rather than the baseline image or the scan as a whole. The ScanProfile tag is only provided on the first, baseline image as it may be large.

Table A.2. TIFF tags

TIFF Tag	Optional?	Description of Content
Software		Starts with “PerkinElmer-QPI”
ImageDescription		Further metadata in XML format (see next section)
ImageWidth		Width of the image in pixels
ImageLength		Height of the image in pixels
ResolutionUnit		Unit used for resolution and position (see below)
XResolution		Pixel X resolution (see below)
YResolution		Pixel Y resolution (see below)

TIFF Tag	Optional?	Description of Content
XPosition	Y	Sample X location in ResolutionUnits. This is ULHC location except for Macro image which reports its image center.
YPosition	Y	Sample Y location in ResolutionUnits. This is ULHC location except for Macro image which reports its image center.
SampleFormat		Integer (1) for BF, FL; or float (3) for extracted components
SMinSampleValue		Minimum signal value in the image
SMaxSampleValue		Maximum signal value in the image
BitsPerSample		8 (FL); 8, 8, 8 (RGB); or 32 (unmixed component)
SamplesPerPixel		1 (FL or unmixed component) or 3 (RGB)
NewSubfileType		0 for full-resolution images, 1 for reduced resolution images
TileWidth	Y	Tile width (512) if tiled format is used
TileLength	Y	Tile height (512) if tiled format is used
TileOffsets	Y	List of tile offsets, if tiled format is used
TileByteCounts	Y	Size of each (compressed) tile, if tiled format is used
StripOffsets	Y	List of strip offsets, if tiled format is not used
RowsPerStrip	Y	Number of rows per strip, if tiled format is not used
StripByteCounts	Y	Size of each (compressed) strip, if tiled format is not used
PlanarConfiguration		1 (chunky) for RGB images, 2 (planar) otherwise.
PlanarConfiguration		1 (chunky) for RGB images, 2 (planar) otherwise
PhotometricInterpretation		2 (RGB) for RGB images, 1 (BlackIsZero) otherwise
DateTime		Acquisition time
Compression		May be None, CCITT Group 3, PackBits, LZW, or JPEG
JPEG fields	Y	JPEG fields are defined when JPEG compression is used

ResolutionUnit, XResolution and YResolution are required fields in a valid TIFF file. When the true resolution of the image is known, ResolutionUnit will be 3 (cm) and XResolution and YResolution will be pixels/cm. When the true resolution is not known, ResolutionUnit will be 2 (inch) and XResolution and YResolution will be 96 (pixels/inch). Pixels from Akoya Biosciences instruments are always square so XResolution and YResolution will always have the same value.

The TIFF spec is not explicit about the data type and value for SMinSampleValue and SMaxSampleValue; the writer uses the same data type as the image pixels (byte or float).

ImageDescription Contents

The ImageDescription tag contains a string in XML format. The string contains a top-level <PerkinElmer-QPI-ImageDescription> element. Nested within this element are child elements with the tag names and values as listed in the table below. Elements appear in the order listed. Values are stored as text content of the element. Elements are required unless otherwise specified. See Table A.3.

Table A.3. ImageDescription Tags

TIFF Tag	Optional?	Contents
DescriptionVersion		Version of the image description field itself, a single number. This document describes version 5 of the field.
AcquisitionSoftware		Software used to acquire the image
Identifier		GUID in string format. This is an identifier for the image file itself.
SlideID	Y	ID of the slide that this image was taken from.
SampleDescription	Y	Text description of the sample
OperatorName	Y	Name of the operator as reported by the operating system running the acquisition software.
Barcode	Y	Barcode text of the slide this image was taken from.
ComputerName	Y	Name of the computer on which the slide was scanned.
ImageType		A string identifying the type of image within the file (Table B.1), with the following values: <ul style="list-style-type: none"> • FullResolution • ReducedResolution • Thumbnail • Overview • Label
IsUnmixedComponent		“True” for unmixed multispectral images, otherwise “False.”
ExposureTime		Exposure time as an integer number of microseconds. For unmixed images, this is the exposure time for the dominant wavelength band for the component (FL); or the brightest wavelength in the cube (BF).
ExposureTimeArray	Y	Exposure times as an integer number of microseconds. Holds multiple values if there is more than one sample per pixel, as for RGB imagery; otherwise the array has a single element.

TIFF Tag	Optional?	Description of contents
Signal Units		<p>A byte <i>www ttt</i> where the <i>ttt</i> nibble indicates the signal unit type from the following:</p> <ul style="list-style-type: none"> – raw counts – normalized (counts/second/gain/full-scale/ binning) – OD (optical density) – dark-corrected counts <p>and the <i>www</i> nibble indicates how the signal is weighted across the spectral bands (or colors):</p> <ul style="list-style-type: none"> – average across all bands – total summed signal across all bands – peak signal in highest-valued band <p>Thus, for example, a value of 68 (hex 44) encodes OD units with peak-signal weighting.</p>
Name	Y	Band name for FL whole slide scans, component name for extracted components. Not present for RGB images
Color	Y	Color to use when rendering this band, as decimal r,g,b byte triplet. Present for FL whole slide scans or extracted components. Not present for RGB images.
Responsivity	Y	Instrument responsivity, if available, for FL whole-slide and unmixed images. See below for details.
Objective	Y	Objective name, if known, otherwise not present. This does not imply a specific pixel resolution – see the XResolution and YResolution standard TIFF tags instead.
ScanProfile	Y	Element containing scan and/or and unmix parameters. It is valid XML whose contents are opaque to most readers. It is only provided on the first (baseline) image and is omitted from all other IFDs.
ProtocolName	Y	Name of the acquisition protocol.
ValidationCode		Used for internal data integrity checks – readers can ignore this.
ImagingCycle	Y	Zero-based acquisition cycle for CODEX imagery
Biomarker	Y	String indicating target such as CD8, Ki67, etc.
ScaleFactor	Y	After allowing for bit-depth conversion, pixel values in this dataset have been scaled by this factor. For ex-ample, if a 12-bit image is stored in an 8-bit image, and counts were scaled from 4095 to 255, this would be 1; if counts were scaled from 3102 to 255, Scale-Factor would be 1.3201. Default is 1.

TIFF Tag	Optional?	Description of contents
NormalizedUnitsFactor	Y	<p>For unmixed data, multiplying the pixel value by this factor will yield component signal in units of normalized counts.</p> <p>Normalized counts are defined by:</p> $Normalized\ counts = \frac{counts}{[FSC * B^2 * G * E * WDF * SW]}$ <p>Where: <i>FSC</i> is the camera full-scale value in counts <i>B</i> is the binning <i>G</i> is the camera gain <i>E</i> is the exposure time, in seconds <i>WDF</i> is a well-depth factor for the sensor used <i>SW</i> is a spectral weighting factor that accounts for the total energy associated with that component's spectrum across all bands.</p>
AutofluorescenceSubtracted	Y	String containing 'true' or 'false', indicating whether this dataset has been subtraction-corrected by a background (autofluorescence) image. Treated as 'false' by default.
LampType	Y	String with the FL lamp type <ul style="list-style-type: none"> - Unknown (default) - PriorLumen200 - PriorLumen220 - XCiteMultiBandLed - PriorMimicLed - SolaLed - Colibri7
BFLampType	Y	String with the BF lamp type <ul style="list-style-type: none"> - Unknown - FlashLamp - WhiteLED - TricolorLED
InstrumentType	Y	String with the instrument type, one of the following: <ul style="list-style-type: none"> - Unknown (default) - Keyence700 - Keyence800 - AxioObserver - Dmi8 - Polaris - Vectra3

TIFF Tag	Optional?	Description of contents
CameraType	Y	String with the camera sensor type, one of the following: <ul style="list-style-type: none"> - SonyICX285AL - HamamatsuC11440 - IMX265 - IMX421
ExcitationFilter	Y	Description of excitation filter with sub-elements listed below If present, all its sub-elements must be provided.
Name		String with the filter public name
Manufacturer		String with the manufacturer or OEM supplier
PartNo		String with the manufacturer part number
Bands		Collection of bands each as described below
Band		Description of a band with child tags noted below
Name		String with the band name
Cuton		Cuton wavelength of this band in nm
Cutoff		Cutoff wavelength of this band in nm
Active		True if this band is excited, False if not
EmissionFilter	Y	Description of emission filter with sub-elements listed below. If present, all its sub-elements must be provided.
Name		String with the filter public name
Manufacturer		String with the manufacturer or EOM supplier
PartNo		String with the manufacturer's part number
Bands		Collection of bands, each described by sub-elements below
Band		Description of a band with sub-elements noted below
Cuton		Cut-on wavelength of this band in nm
Cutoff		Cut-off wavelength of this band in nm
CameraSettings	Y	Description of the camera settings used for this image with sub-elements listed below. If present, all its sub-elements must be provided.
Gain		Gain factor, such as 1, 3, etc.
Binning		Binning extent, such as 1 for full-resolution, 2 for 2x2 binning

TIFF Tag	Optional?	Description of contents
OffsetCounts		Offset count setting for the camera
BitDepth		Bit resolution used to read out the camera. This may differ from the resolution used to store this image (for example, it may be acquired with 12-bit resolution but saved in a QPTIFF file with 8-bit resolution)
Orientation		String containing one of the following values: <ul style="list-style-type: none"> - 'Normal' (default) - 'Rotate180' - 'MirrorV' - 'MirrorH'
ROI	Y	Description of the sensor ROI contributing to the image. Coordinates are sensor pixels before binning. Even if the CameraSettings are provided, this field is optional
X		X origin of ROI. Left-most pixel in sensor array is 0.
Y		Y origin of ROI. Top-most pixel in sensor array is 0.
Width		Width of ROI
Height		Height of ROI

For whole slide images (BF and FL), SignalUnits will be 64 (hex 40) (raw counts, peak signal). For unmixed images, SignalUnits will reflect the unmix settings.

Instrument Responsivity

The <Responsivity> tag is a container for a list of normalized instrument response values. This tag is present for whole-slide FL images from Vectra Polaris and unmixed FL multispectral images originating from Vectra 3 and Vectra Polaris.

For whole-slide images, the <Responsivity> tag will contain one <Filter> tag. The <Filter> tag contents will be different for each image within the TIFF file, reflecting the filter used to take the image.

For unmixed component images, the <Responsivity> tag will contain one <Band> tag for each band in the original image file. The <Band> tags are repeated for each unmixed component image.

The overview, thumbnail and label images do not have <Responsivity> tags.

Each <Filter> or <Band> tag describes the instrument responsivity for acquisitions using that filter or band. Their contents are described in Table A.4.

Table A.4: Contents of <Filter> or <Band> tags within <Responsivity> tag

Tag	Contents
Name	The name of the filter (whole-slide image) or band (component image).
Response	The instrument response to the reference artifact, normalized for exposure. This is raw counts/ ($2^{\text{bit depth}} \times \text{exposure time} \times \text{gain} \times \text{binning area}$), where bit depth is the bit depth of the imagery, exposure time is in seconds, gain is the gain setting of the camera, and binning area is 1 for 1x1, 4 for 2x2, etc.
Date	The date and time of the reference image in UTC, ISO 8601 format.
FilterID	<p>Detailed description of the acquisition filter, as {ExcitationFilter band name}_ {ExcitationFilter manufacturer}:{ExcitationFilter part number} / {EmissionFilter manufacturer}:{EmissionFilter part number}</p> <p>The {ExcitationFilter band name}_ prefix may be omitted if there is only a single excitation band.</p>

Sample ImageDescription

Sample ImageDescription for the DAPI band of a FL whole-slide image, containing a single <Filter> tag:

```
<?xml version="1.0" encoding="utf-8"?>
< PerkinElmer-QPI-ImageDescription >
  <DescriptionVersion>1</DescriptionVersion>
  <AcquisitionSoftware>VectraScan 1.0.0</AcquisitionSoftware>
  <ImageType>FullResolution</ImageType>
  <Identifier>AABED946-BB58-44FB-95B3-48E177E3BB83</Identifier>
  <IsUnmixedComponent>False</IsUnmixedComponent>
  <ExposureTime>50</ExposureTime>
  <SignalUnits>64</SignalUnits>
  <Name>DAPI</Name>
  <Color>0,0,255</Color>
  <Responsivity>
    <Filter>
      <Name>DAPI</Name>
      <Response>30.7</Response>
      <Date>2015-10-22T13:10:18.0618849Z</Date>
      <FilterID>Semrock:FF02-409/LP-25 Emission / Semrock:FF01-377/50-25
        Excitation</FilterID>
    </Filter>
  </Responsivity>
  <Objective>4x</Objective>
  <ScanProfile><!-- this will be a serialized scan protocol. It is valid XML but otherwise
    opaque -->
  </ScanProfile>
  <ValidationCode>4281ff86778db65892c05151d5de738d</ValidationCode>
</ PerkinElmer-QPI-ImageDescription >
```

Appendix B: Akoya Software & Spectra Compatibility

Some files are read by multiple Akoya acquisition and analysis software programs to facilitate unmixing and analysis workflows. Due to major updates in the combined release of Phenolmager HT 2.0.0, Phenochart 2.0.0, and inForm 3.0.0 software, outputs from these versions are not backward compatible with previous versions of software. However, files from previous versions of software are forward compatible. Software compatibility considerations are detailed in Table B.1.

Table B.1. Akoya Software Compatibility Considerations

	VectraPolaris 1.x.x	Phenolmager HT 2.x.x	Phenochart 1.x.x	Phenochart 2.x.x	inForm 2.x.x	inForm 3.x.x
Scanning protocols created in VectraPolaris 1.x.x can be opened by	Yes	Yes	N/A	N/A	N/A	N/A
Scanning protocols created in Phenolmager HT 2.x.x can be opened by	No	Yes	N/A	N/A	N/A	N/A
Imagery acquired with VectraPolaris 1.x.x can be opened by	N/A	N/A	Yes	Yes	Yes	Yes
Imagery acquired with Phenolmager HT 2.x.x can be opened by	N/A	N/A	Some*	Yes	No	Yes
Annotations created by Phenochart 1.x.x can be opened by	Yes	Yes	Yes	Yes	Yes	Yes
Annotations created by Phenochart 2.x.x can be opened by	No	Yes	No	Yes	No	Yes
Algorithms created by inForm 2.x.x can be opened by	N/A	N/A	N/A	N/A	Yes	Yes
Algorithms created by inForm 3.x.x can be opened by	N/A	N/A	N/A	N/A	No	Yes
Projects created by inForm 2.x.x can be opened by	N/A	N/A	N/A	N/A	Yes	Yes**
Projects created by inForm 3.x.x can be opened by	N/A	N/A	N/A	N/A	No	Yes

	VectraPolaris 1.x.x	PhenolMager HT 2.x.x	Phenochart 1.x.x	Phenochart 2.x.x	inForm 2.x.x	inForm 3.x.x
Spectra extracted in inForm 2.x.x can be exported as an un-mixing library for HT 2.x.x by	N/A	N/A	N/A	N/A	No	Yes
Spectra extracted in inForm 3.x.x can be exported as an un-mixing library for HT 2.x.x by	N/A	N/A	N/A	N/A	No	Yes

*Can be opened, but will not allow for stamping for inForm.

**Background on Projects created by inForm 2.x opening in inForm 3.x. Two considerations:

1. Beginning with inForm 3.0, there are slight changes to cell splitting or segmentation due to the different numerical routines used. Also, the cell numbering can be assigned differently even for the same segmentation results. These do not affect the results in a significant way for normal use, but as a result, the assignment of cell ID #s to cells can differ between 2.x and 3.x. This means that phenotyping training sets from inForm 2.x projects cannot be used for training phenotyping classifiers in 3.x because the training sets consist of tables of phenotypes keyed by cell ID. If the cell ID is wrong, the training data is erroneous so the resulting classifier will not be properly trained. For this reason, inForm 3.x will not bring in training cell ID information for projects from inForm 2.x. It will post this message when this happens: "[Project] was created with inForm 2.8 or earlier, and phenotype training data from 2.8 cannot be imported into inForm 3.0. All cell phenotype training data has been discarded for that reason. The algorithm's phenotype step can be run with its current training, or use inForm 2.8 or earlier to add training cells or revise phenotyping."
2. For projects created in inForm 2.8 with imagery that has 10+ channels (i.e., PhenoCycler imagery), phenotyping computation may become unusably slow in inForm 3.x. It is recommended that inForm 2.8 projects with imagery that has 10+ channels are only run in inForm 2.8. New projects created in inForm 3.0 do not have this issue.

Considerations for Side-by-Side Installations

- Vectra Polaris 1.x and PhenolMager HT 2.x software can be installed side-by-side. If needed, Vectra Polaris 1.x may continue to be used for some ongoing studies that use inForm 2.x algorithms.
- Discerning which version is launching will be easy due to differing program names. Phenochart 1.x and 2.x can be installed side-by-side.
- Users often launch Phenochart by double-clicking a QPTIFF. In this case, the launched version might not be obvious but will be easily confirmed via About Box. inForm 2.x and 3.x can be installed side-by-side.
- Users sometimes launch inForm by double-clicking on a project. In this case, the launched version might not be obvious but will be easily confirmed via About Box.

Spectra compatibility considerations are detailed in Table B.2.

Table B.2. Spectra compatibility and incompatibility considerations

	Unmix 8-bit QTIFFs in inForm	Unmix Extended-bit QTIFFs in inForm	Unmix in HT via .umxlib
Synthetic Opal & AP spectra can be applied to	YES	YES	YES
Opal spectra extracted from 8-bit QTIFFs can be applied to	YES	YES	YES
AF spectra extracted from 8-bit QTIFFs can be applied to	YES	NO	YES
Opal spectra extracted from Extended-bit QTIFFs can be applied to	YES	YES	YES
AF spectra extracted from Extended-bit QTIFFs can be applied to	NO	YES	YES

Appendix C: Phenolmager HT EULA

The following is an agreement (the "Agreement") between you and Akoya Biosciences, Inc., 100 Campus Drive, 6th Floor, Marlborough, MA 01752 ("Akoya") for software known as Research Use Only Phenolmager HT, including any updates, upgrades or enhancements provided by Akoya and any and all accompanying documentation (collectively, the "Software").

AKOYA IS WILLING TO LICENSE THE SOFTWARE TO YOU ONLY UPON THE CONDITION THAT YOU ACCEPT ALL OF THE TERMS CONTAINED IN THIS AGREEMENT. PLEASE READ THE TERMS CAREFULLY BEFORE USING THE SOFTWARE, AS INSTALLING AND/OR USING THE SOFTWARE WILL INDICATE YOUR AGREEMENT TO THEM. IF YOU DO NOT AGREE TO ALL OF THE TERMS AND CONDITIONS IN THIS AGREEMENT, AKOYA IS NOT WILLING TO LICENSE THE SOFTWARE TO YOU, AND YOU (1) MAY NOT USE THE SOFTWARE AND (2) MUST CONTACT AKOYA WITHIN THREE (3) DAYS OF RECEIPT TO ARRANGE FOR THE RETURN OF THE SOFTWARE. BY INSTALLING AND/OR USING THE SOFTWARE, YOU ALSO REPRESENT AND WARRANT THAT YOU HAVE THE AUTHORITY TO AGREE TO AND BE BOUND BY THE TERMS OF THE AGREEMENT.

1. Single Use License. The Software is licensed to you and not sold, and you acknowledge that there is no transfer of title or ownership of the Software to you. Subject to the terms and conditions of this Agreement and your payment of the license fees for the Software, Akoya hereby grants to you a restricted, non-exclusive, non-transferable, non-assignable, non-sublicensable and revocable single-use license to use, for your internal research purposes only, the executable code version of the Software and the accompanying documentation in hard copy or electronic format. Akoya grants no other licenses hereunder, except for the Software set forth above. Nothing in this Agreement shall be construed as a general internal use license. For the avoidance of doubt, no right to use the Software in diagnostic applications is granted under this Agreement. AKOYA RESERVES ALL RIGHTS NOT EXPRESSLY GRANTED BY THIS AGREEMENT.

2. Specific Restrictions. You may use the Software only on a single computer. You may make only one (1) copy of the Software, solely for backup purposes, and such copy must reproduce and include all copyright notices indicating Akoya's ownership of the Software. Except as expressly permitted by applicable law or Akoya, you agree that neither you nor a third party acting on your behalf, directly or indirectly, shall: (i) decompile, disassemble or reverse engineer the Software or otherwise attempt to derive the source code of the Software or underlying ideas, techniques or algorithms, including the review of data structures or similar materials produced by the Software; (ii) modify, alter, translate, or create derivative works of the Software; (iii) lease, rent, loan, commercially share or otherwise use the Software in any manner to provide service bureau, commercial time-sharing or other computer services to third parties; (iv) transmit the Software or provide its functionality, in whole or in part, over the Internet or other network; (v) use the Software for any purpose other than solely for your internal research; (vi) use the Software to provide any analytics or diagnostics or otherwise for the benefit of any third party; (vii) sell, distribute, publish, rent, lease, sublicense, assign possession, allow use of or otherwise transfer the Software or any copy of the Software to a third party; (viii) create any computer software program or documentation which is similar to or competitive with the Software; (ix) encumber, disclose or otherwise make available the Software to any third party; or (x) remove, modify or obscure any copyright, trademark, or other proprietary rights notice in or on the Software. This Agreement may not be assigned in whole or in part, except that either party may assign this Agreement without the other party's prior written consent in the event of a merger with an affiliate or other

reorganization with an affiliate. Any attempted assignment contrary to the terms of this Agreement shall be void.

3. Ownership of Software. As between the parties, Akoya owns all right, title and interest, including all existing or future copyrights, trademarks, tools, know-how, and processes, in and to the Software and any copies, modifications, alterations or derivative works thereof. You shall have no title to or ownership of the Software. As between the parties, Akoya shall retain all rights to market, license and/or sell the Software. The Software (including accompanying documentation) constitutes the unpublished confidential information and trade secrets of Akoya. You (i) acknowledge, and agree not to contest, Akoya's rights in the Software; and (ii) agree not to disclose and to take reasonable steps to protect any and all confidential information of Akoya regarding the Software or that is otherwise disclosed to you in connection with this Agreement. You agree that you are responsible for ensuring your employees and agents are prohibited from taking any action otherwise prohibited to you under this Agreement and to notify Akoya promptly and in writing of the circumstances surrounding any possession, use or knowledge of the Software or any part thereof known to you which is prohibited hereunder. You recognize that the covenants contained in this Agreement are reasonable and necessary to protect the legitimate interests of Akoya, that Akoya would not have entered into this Agreement in the absence of such covenants, and that your breach or threatened breach of such covenants shall cause Akoya irreparable harm and significant injury, the amount of which shall be extremely difficult to estimate and ascertain, thus making any remedy at law or in damages inadequate. Therefore, you agree that Akoya shall be entitled, without the necessity of posting any bond or security, to the issuance of injunctive relief by any court of competent jurisdiction enjoining any breach or threatened breach of such covenants and for any other relief such court deems appropriate. This right shall be in addition to any other remedy available to Akoya at law or in equity.

4. Termination. The license granted in Section 1 above is effective until terminated. This Agreement is conditioned upon your continued compliance with the terms and conditions hereof and will terminate automatically without notice from Akoya if you fail to comply with any term or condition of this Agreement. Furthermore, Akoya may terminate this Agreement at any time upon thirty (30) days' notice. Upon termination of this Agreement, your license to use the Software or any part thereof shall immediately terminate, and you shall immediately destroy all copies of the Software (including all accompanying documentation) and any other confidential and proprietary information you have received during or in connection with this Agreement. Termination of this Agreement by Akoya shall not entitle you to a refund of any fees paid for the license of the Software.

5. Limited Warranty. The sole and exclusive warranty applicable to software and firmware products provided by Akoya for use with a processor is as follows: Akoya warrants that such software and firmware will substantially conform to Akoya's documentation and program manuals current at the time of shipment to you when properly installed on the processor. The limited warranty set forth above shall not apply if the claimed problem or issue is caused by: the malfunction of computer hardware or software not provided by Akoya; your negligence or fault; your failure to follow the instructions set forth in the documentation; changes in the hardware or other software not authorized by Akoya; modifications or changes to the Software not made or suggested by Akoya; or your failure to implement and maintain a proper backup and recovery system. Your exclusive remedy under this Section 5 shall be, at Akoya's option, the replacement or correction of the problem or issue causing the non-conformance with the above warranty as promptly as possible following notice of such non-conformity, or a refund of the price paid for the Software.

6. No Other Warranties. EXCEPT FOR THE LIMITED WARRANTY STATED IMMEDIATELY ABOVE AND TO THE MAXIMUM EXTENT PERMITTED BY APPLICABLE LAW, THE SOFTWARE IS PROVIDED "AS IS" WITHOUT WARRANTY OF ANY KIND, AND AKOYA EXPRESSLY DISCLAIMS ANY AND ALL IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY WARRANTY OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NON-INFRINGEMENT OF THIRD PARTY RIGHTS. AKOYA DOES NOT REPRESENT OR WARRANT THAT THE RESULTS OR THE USE OF THE SOFTWARE WILL BE CORRECT, ACCURATE OR RELIABLE, WILL MEET YOUR REQUIREMENTS, OR THAT THE SOFTWARE WILL OPERATE UNINTERRUPTED OR ERROR-FREE, OR THAT DEFECTS IN THE SOFTWARE WILL BE CORRECTED. Some jurisdictions do not allow the exclusion of implied warranties, so the above disclaimer may not apply to you. YOU ASSUME ALL RISK ASSOCIATED WITH THE USE, RESULTS, AND PERFORMANCE OF THE SOFTWARE.

7. Limitation of Liability. IN NO EVENT SHALL AKOYA, ITS AFFILIATES OR SUPPLIERS, OR THEIR RESPECTIVE EMPLOYEES, OFFICERS OR AGENTS, BE LIABLE FOR ANY DAMAGES ARISING OUT OF THE DELIVERY, INSTALLATION, PERFORMANCE, OR USE OF OR INABILITY TO USE THE SOFTWARE, INCLUDING WITHOUT LIMITATION INCIDENTAL, SPECIAL, CONSEQUENTIAL, PUNITIVE, EXEMPLARY, INDIRECT OR DIRECT DAMAGES (INCLUDING WITHOUT LIMITATION DAMAGES FOR LOSS OF PROFITS, LOSS OF DATA, RE-RUN TIME, INACCURATE INPUT, WORK DELAYS, BUSINESS INTERRUPTION OR ANY OTHER COMMERCIAL DAMAGES OR LOSSES), WHETHER IN AN ACTION IN CONTRACT, TORT (INCLUDING NEGLIGENCE AND STRICT LIABILITY) OR OTHERWISE, AND EVEN IF AKOYA HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES, AND REGARDLESS OF WHETHER ANY REMEDY FAILS OF ITS ESSENTIAL PURPOSE. AKOYA'S AGGREGATE LIABILITY ARISING OUT OF OR RELATING TO THIS AGREEMENT, FOR ANY AND ALL DAMAGES THAT YOU MIGHT INCUR (REGARDLESS OF THE FORM OF ACTION GIVING RISE TO SUCH LIABILITY, WHETHER IN CONTRACT, TORT OR OTHERWISE) SHALL NOT EXCEED THE PRICE PAID BY YOU FOR THE LICENSE OF THE SOFTWARE. Some jurisdictions do not allow the exclusion or limitation of liability for consequential or incidental damages. In such jurisdictions, our liability is limited to the greatest extent permitted by law.

8. U.S. Government End Users. The Software qualifies as commercial computer software for purposes of the Federal Acquisition Regulations (FAR) §52.227-19 and the Department of Defense Supplement to the FAR (DFARS) §52.227-7013. If the Software is acquired by a civilian government agency, it is furnished with only the minimum Restricted Rights provided by FAR 52.227-19. If the Software is acquired by a military agency, it is furnished with only the minimum Restricted Rights provided by DFARS 52.227-7013(c)(1)(ii).

9. Miscellaneous. This Agreement contains the entire agreement of the parties with respect to the subject matter hereof and supersedes any proposal or prior agreement, written or oral, and any other communications between the parties relating to the subject matter hereof. No modification or waiver of any provision of this Agreement shall be effective unless in writing and signed by the parties. No delay or failure on the part of any party in exercising any right under this Agreement shall impair any such right or any remedy of such party, nor shall it be construed to be a waiver of any continuing breach or default under this Agreement. No action, regardless of form, arising out of the transactions contemplated herein may be brought by either party more than one (1) year after the cause of action has accrued. In the event any provision of this Agreement is held to be unenforceable, the remaining provisions of this Agreement will remain in full force and effect. This Agreement shall be governed by the laws of the State of Delaware without regard to principles of conflicts of laws. Any disputes relating hereto shall be adjudicated only in any court of competent jurisdiction in the State of Delaware, and you hereby consent to the exclusive jurisdiction of those courts. This Agreement shall not be governed by the United Nations Convention on Contracts for the International Sale of Goods, the application of which is

expressly excluded. You may not assign or otherwise transfer this Agreement or any of your rights or obligations therein without the prior written consent of Akoya. You may not use the Software for any unlawful purpose nor export or re-export the Software except as authorized by law.

10. This software uses the NLog 5.2.5 libraries (<https://github.com/NLog/NLog>)

This software makes use of the NLog libraries. These libraries have not been modified by Akoya and all rights are reserved by the copyright holder. This software is provided "as is" without express or implied warranty, and with no claim as to its suitability for any purpose. The libraries are covered by the BSD license.

The BSD for NLog license is reproduced below:

Copyright (c) 2004-2021 Jaroslaw Kowalski <jaak@jkwowski.net>, Kim Christensen, Julian Verdurmen. All rights reserved.

Redistribution and use in source and binary forms, with or without modification, are permitted provided that the following conditions are met:

* Redistributions of source code must retain the above copyright notice, this list of conditions and the following disclaimer.

* Redistributions in binary form must reproduce the above copyright notice, this list of conditions and the following disclaimer in the documentation and/or other materials provided with the distribution.

* Neither the name of Jaroslaw Kowalski nor the names of its contributors may be used to endorse or promote products derived from this software without specific prior written permission.

THIS SOFTWARE IS PROVIDED BY THE COPYRIGHT HOLDERS AND CONTRIBUTORS "AS IS" AND ANY EXPRESS OR IMPLIED WARRANTIES, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE ARE DISCLAIMED. IN NO EVENT SHALL THE COPYRIGHT OWNER OR CONTRIBUTORS BE LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, EXEMPLARY, OR CONSEQUENTIAL DAMAGES (INCLUDING, BUT NOT LIMITED TO, PROCUREMENT OF SUBSTITUTE GOODS OR SERVICES; LOSS OF USE, DATA, OR PROFITS; OR BUSINESS INTERRUPTION) HOWEVER CAUSED AND ON ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, STRICT LIABILITY, OR TORT (INCLUDING NEGLIGENCE OR OTHERWISE) ARISING IN ANY WAY OUT OF THE USE OF THIS SOFTWARE, EVEN IF ADVISED OF THE POSSIBILITY OF SUCH DAMAGE.

11. This software uses the LibTIFF 4.0.6 library (<http://remotesensing.org/libtiff>)

Copyright (c) 1988-1997 Sam Leffler Copyright (c) 1991-1997 Silicon Graphics, Inc.

Permission to use, copy, modify, distribute, and sell this software and its documentation for any purpose is hereby granted without fee, provided that (i) the above copyright notices and this permission notice appear in all copies of the software and related documentation, and (ii) the names of Sam Leffler and Silicon Graphics may not be used in any advertising or publicity relating to the software without the specific, prior written permission of Sam Leffler and Silicon Graphics.

THE SOFTWARE IS PROVIDED "AS-IS" AND WITHOUT WARRANTY OF ANY KIND, EXPRESS, IMPLIED OR OTHERWISE, INCLUDING WITHOUT LIMITATION, ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

IN NO EVENT SHALL SAM LEFFLER OR SILICON GRAPHICS BE LIABLE FOR ANY SPECIAL, INCIDENTAL, INDIRECT OR CONSEQUENTIAL DAMAGES OF ANY KIND, OR ANY DAMAGES WHATSOEVER RESULTING FROM LOSS OF USE, DATA OR PROFITS, WHETHER OR NOT ADVISED OF THE POSSIBILITY OF DAMAGE, AND ON ANY THEORY OF LIABILITY, ARISING OUT OF OR IN CONNECTION WITH THE USE OR PERFORMANCE OF THIS SOFTWARE.

12. This software uses the Vigna 1.6.0 library (<http://hci.iwr.uni-heidelberg.de/vigna>)

This software includes machine-executable object code generated by a source language processor from the Vigna libraries covered by the VIGRA license. These libraries have not been modified by Akoya and all rights are reserved by the copyright holder. This software is provided "as is" without express or implied warranty, and with no claim as to its suitability for any purpose.

The Vigna library is Copyright © 2008, Ulrich Kothe, Heidelberg Collaboratory for Image Processing

The VIGRA License

=====

(identical to the MIT X11 License)

Permission is hereby granted, free of charge, to any person obtaining a copy of this software and associated documentation files (the "Software"), to deal in the Software without restriction, including without limitation the rights to use, copy, modify, merge, publish, distribute, sublicense, and/or sell copies of the Software, and to permit persons to whom the Software is furnished to do so, subject to the following conditions:

The above copyright notice and this permission notice shall be included in all copies or substantial portions of the Software.

THE SOFTWARE IS PROVIDED "AS IS", WITHOUT WARRANTY OF ANY KIND EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND NONINFRINGEMENT. IN NO EVENT SHALL THE AUTHORS OR COPYRIGHT HOLDERS BE LIABLE FOR ANY CLAIM, DAMAGES OR OTHER LIABILITY, WHETHER IN AN ACTION OF CONTRACT, TORT OR OTHERWISE, ARISING FROM, OUT OF OR IN CONNECTION WITH THE SOFTWARE OR THE USE OR OTHER DEALINGS IN THE SOFTWARE.

13. This software uses the SMIL 1.0.0 library (<https://smil.cmm.minesparis.psl.eu/doc>)

This software makes use of the Autofac library. This library has not been modified by Akoya and all rights are reserved by the copyright holder. This software is provided "as is" without express or implied warranty, and with no claim as to its suitability for any purpose. The library is covered by the license reproduced below:

Copyright (c) 2011-2016 Matthieu FAESSEL and ARMINES

Copyright (c) 2017-2022 CMM - Centre de Morphologie Mathematique -Mines-ParisTech

All rights reserved.

Redistribution and use in source and binary forms, with or without modification, are permitted provided that the following conditions are met:

- Redistributions of source code must retain the above copyright notice, this list of conditions and the following disclaimer.
- Redistributions in binary form must reproduce the above copyright notice, this list of conditions and the following disclaimer in the documentation and/or other materials provided with the distribution.
- Neither the name of Matthieu FAESSEL, or ARMINES nor the names of its contributors may be used to endorse or promote products derived from this software without specific prior written permission.

THIS SOFTWARE IS PROVIDED BY THE COPYRIGHT HOLDERS AND CONTRIBUTORS "AS IS" AND ANY EXPRESS OR IMPLIED WARRANTIES, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE

ARE DISCLAIMED. IN NO EVENT SHALL THE COPYRIGHT HOLDERS AND CONTRIBUTORS BE LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, EXEMPLARY, OR CONSEQUENTIAL DAMAGES (INCLUDING, BUT NOT LIMITED TO, PROCUREMENT OF SUBSTITUTE GOODS OR SERVICES; LOSS OF USE, DATA, OR PROFITS; OR BUSINESS INTERRUPTION) HOWEVER CAUSED AND ON ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, STRICT LIABILITY, OR TORT (INCLUDING NEGLIGENCE OR OTHERWISE) ARISING IN ANY WAY OUT OF THE USE OF THIS SOFTWARE, EVEN IF ADVISED OF THE POSSIBILITY OF SUCH DAMAGE.

14. This software uses the libjpeg-turbo 1.5.0 library (<http://libjpeg-turbo.virtualgl.org>)

This software makes use of the libjpeg-turbo library. This software is based in part on the work of the Independent JPEG Group. This library has not been modified by Akoya and all rights are reserved by the copyright holder. This software is provided "as is" without express or implied warranty, and with no claim as to its suitability for any purpose. The library is covered by IJG License. Additionally libjpeg-turbo incorporates parts of libjpeg which is covered by the BSD license.

The IJG license is reproduced below:

The authors make NO WARRANTY or representation, either express or implied, with respect to this software, its quality, accuracy, merchantability, or fitness for a particular purpose. This software is provided "AS IS", and you, its user, assume the entire risk as to its quality and accuracy.

This software is copyright (C) 1991-1998, Thomas G. Lane. All Rights Reserved except as specified below.

Permission is hereby granted to use, copy, modify, and distribute this software (or portions thereof) for any purpose, without fee, subject to these conditions:

(1) If any part of the source code for this software is distributed, then this README file must be included, with this copyright and no-warranty notice unaltered; and any additions, deletions, or changes to the original files must be clearly indicated in accompanying documentation.

(2) If only executable code is distributed, then the accompanying documentation must state that "this software is based in part on the work of the Independent JPEG Group".

(3) Permission for use of this software is granted only if the user accepts full responsibility for any undesirable consequences; the authors accept NO LIABILITY for damages of any kind. These conditions apply to any software derived from or based on the IJG code, not just to the unmodified library. If you use our work, you ought to acknowledge us.

Permission is NOT granted for the use of any IJG author's name or company name in advertising or publicity relating to this software or products derived from it. This software may be referred to only as "the Independent JPEG Group's software".

We specifically permit and encourage the use of this software as the basis of commercial products, provided that all warranty or liability claims are assumed by the product vendor. ansi2knr.c is included in this distribution by permission of L. Peter Deutsch, sole proprietor of its copyright holder, Aladdin Enterprises of Menlo Park, CA. ansi2knr.c is NOT covered by the above copyright and conditions, but instead by the usual distribution terms of the Free Software Foundation; principally, that you must include source code if you redistribute it. (See the file ansi2knr.c for full details.) However, since ansi2knr.c is not needed as part of any program generated from the IJG code, this does not limit you more than the foregoing paragraphs do.

The Unix configuration script "configure" was produced with GNU Autoconf. It is copyright by the Free Software Foundation but is freely distributable. The same holds for its supporting scripts (config.guess, config.sub, ltconfig, ltmain.sh). Another support script,

install-sh, is copyright by M.I.T. but is also freely distributable.

It appears that the arithmetic coding option of the JPEG spec is covered by patents owned by IBM, AT&T, and Mitsubishi. Hence arithmetic coding cannot legally be used without obtaining one or more licenses. For this reason, support for arithmetic coding has been removed from the free JPEG software. (Since arithmetic coding provides only a marginal gain over the unpatented Huffman mode, it is unlikely that very many implementations will support it.) So far as we are aware, there are no patent restrictions on the remaining code.

The IJG distribution formerly included code to read and write GIF files. To avoid entanglement with the Unisys LZW patent, GIF reading support has been removed altogether, and the GIF writer has been simplified to produce "uncompressed GIFs". This technique does not use the LZW algorithm; the resulting GIF files are larger than usual, but are readable by all standard GIF decoders.

We are required to state that "The Graphics Interchange Format(c) is the Copyright property of CompuServe Incorporated. GIF(sm) is a Service Mark property of CompuServe Incorporated."

The BSD for libjpeg license is reproduced below:

The authors make NO WARRANTY or representation, either express or implied, with respect to this software, its quality, accuracy, merchantability, or fitness for a particular purpose. This software is provided "AS IS", and you, its user, assume the entire risk as to its quality and accuracy.

This software is copyright (C) 1991-2010, Thomas G. Lane, Guido Vollbeding.

All Rights Reserved except as specified below.

Permission is hereby granted to use, copy, modify, and distribute this software (or portions thereof) for any purpose, without fee, subject to these conditions:

(1) If any part of the source code for this software is distributed, then this README file must be included, with this copyright and no-warranty notice unaltered; and any additions, deletions, or changes to the original files must be clearly indicated in accompanying documentation.

(2) If only executable code is distributed, then the accompanying documentation must state that "this software is based in part on the work of the Independent JPEG Group".

(3) Permission for use of this software is granted only if the user accepts full responsibility for any undesirable consequences; the authors accept NO LIABILITY for damages of any kind. These conditions apply to any software derived from or based on the IJG code, not just to the unmodified library. If you use our work, you ought to acknowledge us.

Permission is NOT granted for the use of any IJG author's name or company name in advertising or publicity relating to this software or products derived from it. This software may be referred to only as "the Independent JPEG Group's software".

We specifically permit and encourage the use of this software as the basis of commercial products, provided that all warranty or liability claims are assumed by the product vendor.

15. This software uses the zlib 1.2.8 library (<http://www.zlib.net>)

This software makes use of the zlib library. This library has not been modified by Akoya and all rights are reserved by the copyright holder. This software is provided "as is" without express or implied warranty, and with no claim as to its suitability for any purpose. The library is covered by the zlib License.

The zlib 1 license is reproduced below:

(C) 1995-2012 Jean-loup Gailly and Mark Adler

This software is provided 'as-is', without any express or implied warranty. In no event will the authors be held liable for any damages arising from the use of this software.

Permission is granted to anyone to use this software for any purpose, including commercial applications, and to alter it and redistribute it freely, subject to the following restrictions:

1. The origin of this software must not be misrepresented; you must not claim that you wrote the original software. If you use this software in a product, an acknowledgment in the product documentation would be appreciated but is not required.
2. Altered source versions must be plainly marked as such, and must not be misrepresented as being the original software.
3. This notice may not be removed or altered from any source distribution.

16. This software uses FTDI drivers (<http://www.ftdichip.com>)

This software makes use of drivers provided by Future Technology Devices International Limited, which are "as is" and any express or implied warranties, including, but not limited to, the implied warranties of merchantability and fitness for a particular purpose are disclaimed. In no event shall future technology devices international limited be liable for any direct, indirect, incidental, special, exemplary, or consequential damages (including, but not limited to, procurement of substitute goods or services; loss of use, data, or profits; or business interruption) however caused and on any theory of liability, whether in contract, strict liability, or tort (including negligence or otherwise) arising in any way out of the use of this software, even if advised of the possibility of such damage. FTDI drivers may be used only in conjunction with products based on FTDI parts. FTDI drivers may be distributed in any form as long as license information is not modified.

Copyright Future Technology Devices International Limited of 2 Seaward Place, Centurion Business Park, Glasgow, Scotland, G41 1HH

17. This software uses the xzing.NET 0.16.9 libraries

(<https://www.nuget.org/packages/ZXing.Net>)

This software makes use of the xzing.NET libraries. These libraries have not been modified by Akoya and all rights are reserved by the copyright holder. This software is provided "as is" without express or implied warranty, and with no claim as to its suitability for any purpose.

Copyright 2022 Michael Jahn

Licensed under the Apache License, Version 2.0 (the "License"); you may not use this file except in compliance with the License. You may obtain a copy of the License at <http://www.apache.org/licenses/LICENSE-2.0>

Unless required by applicable law or agreed to in writing, software distributed under the License is distributed on an "AS IS" BASIS, WITHOUT WARRANTIES OR CONDITIONS OF ANY KIND, either express or implied. See the License for the specific language governing permissions and limitations under the License.

18. This software uses the Boost 1.68.0 library (<http://www.boost.org>)

This software includes machine-executable object code generated by a source language processor from the Boost libraries covered by the Boost license. These libraries have not been modified by Akoya and all rights are reserved by the copyright holder. This software is provided "as is" without express or implied warranty, and with no claim as to its suitability for any purpose.

The Boost license is reproduced below:

Boost Software License - Version 1.0 - August 17th, 2003

Permission is hereby granted, free of charge, to any person or organization obtaining a copy of the software and accompanying documentation covered by this license (the "Software") to use, reproduce, display, distribute, execute, and transmit the Software, and to prepare derivative works of the Software, and to permit third-parties to whom the Software is furnished to do so, all subject to the following:

The copyright notices in the Software and this entire statement, including the above license grant, this restriction and the following disclaimer, must be included in all copies of the Software, in whole or in part, and all derivative works of the Software, unless such copies or derivative works are solely in the form of machine-executable object code generated by a source language processor.

THE SOFTWARE IS PROVIDED "AS IS", WITHOUT WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, TITLE AND NON-INFRINGEMENT. IN NO EVENT SHALL THE COPYRIGHT HOLDERS OR ANYONE DISTRIBUTING THE SOFTWARE BE LIABLE FOR ANY DAMAGES OR OTHER LIABILITY, WHETHER IN CONTRACT, TORT OR OTHERWISE, ARISING FROM, OUT OF OR IN CONNECTION WITH THE SOFTWARE OR THE USE OR OTHER DEALINGS IN THE SOFTWARE.

19. This software uses the Fugue icons, version 3.5.6 (<http://p.yusukekamiyamane.com>)

This software makes use of the Fugue icons. This content has not been modified by Akoya and all rights are reserved by the copyright holder. This software is provided "as is" without express or implied warranty, and with no claim as to its suitability for any purpose. The icons are covered by the Creative Commons 3.0 License.

(C) 2013 Yusuke Kamiyamane. All rights reserved.

The Creative Commons 3.0 License is reproduced below:

THE WORK (AS DEFINED BELOW) IS PROVIDED UNDER THE TERMS OF THIS CREATIVE COMMONS PUBLIC LICENSE ("CCPL" OR "LICENSE"). THE WORK IS PROTECTED BY COPYRIGHT AND/OR OTHER APPLICABLE LAW. ANY USE OF THE WORK OTHER THAN AS AUTHORIZED UNDER THIS LICENSE OR COPYRIGHT LAW IS PROHIBITED.

BY EXERCISING ANY RIGHTS TO THE WORK PROVIDED HERE, YOU ACCEPT AND AGREE TO BE BOUND BY THE TERMS OF THIS LICENSE. TO THE EXTENT THIS LICENSE MAY BE CONSIDERED TO BE A CONTRACT, THE LICENSOR GRANTS YOU THE RIGHTS CONTAINED HERE IN CONSIDERATION OF YOUR ACCEPTANCE OF SUCH TERMS AND CONDITIONS.

1. Definitions

- a. **"Adaptation"** means a work based upon the Work, or upon the Work and other pre-existing works, such as a translation, adaptation, derivative work, arrangement of music or other alterations of a literary or artistic work, or phonogram or performance and includes cinematographic adaptations or any other form in which the Work may be recast, transformed, or adapted including in any form recognizably derived from the original, except that a work that constitutes a Collection will not be considered an Adaptation for the purpose of this License. For the avoidance of doubt, where the Work is a musical work, performance or phonogram, the synchronization of the Work in timed-relation with a moving image ("synching") will be considered an Adaptation for the purpose of this License.
- b. **"Collection"** means a collection of literary or artistic works, such as encyclopedias and anthologies, or performances, phonograms or broadcasts, or other works or subject matter other than works listed in Section 1(f) below, which, by reason of the selection and arrangement of their contents, constitute intellectual creations, in which the Work is included in its entirety in unmodified form along with one or more other contributions, each constituting separate and independent works in themselves, which together are assembled into a collective whole. A work that constitutes a Collection will not be considered an Adaptation (as defined above) for the purposes of this License.
- c. **"Distribute"** means to make available to the public the original and copies of the Work or Adaptation, as appropriate, through sale or other transfer of ownership.
- d. **"Licensor"** means the individual, individuals, entity or entities that offer(s) the

- Work under the terms of this License.
- e. **"Original Author"** means, in the case of a literary or artistic work, the individual, individuals, entity or entities who created the Work or if no individual or entity can be identified, the publisher; and in addition (i) in the case of a performance the actors, singers, musicians, dancers, and other persons who act, sing, deliver, declaim, play in, interpret or otherwise perform literary or artistic works or expressions of folklore; (ii) in the case of a phonogram the producer being the person or legal entity who first fixes the sounds of a performance or other sounds; and, (iii) in the case of broadcasts, the organization that transmits the broadcast.
 - f. **"Work"** means the literary and/or artistic work offered under the terms of this License including without limitation any production in the literary, scientific and artistic domain, whatever may be the mode or form of its expression including digital form, such as a book, pamphlet and other writing; a lecture, address, sermon or other work of the same nature; a dramatic or dramatico-musical work; a choreographic work or entertainment in dumb show; a musical composition with or without words; a cinematographic work to which are assimilated works expressed by a process analogous to cinematography; a work of drawing, painting, architecture, sculpture, engraving or lithography; a photographic work to which are assimilated works expressed by a process analogous to photography; a work of applied art; an illustration, map, plan, sketch or three-dimensional work relative to geography, topography, architecture or science; a performance; a broadcast; a phonogram; a compilation of data to the extent it is protected as a copyrightable work; or a work performed by a variety or circus performer to the extent it is not otherwise considered a literary or artistic work.
 - g. **"You"** means an individual or entity exercising rights under this License who has not previously violated the terms of this License with respect to the Work, or who has received express permission from the Licensor to exercise rights under this License despite a previous violation.
 - h. **"Publicly Perform"** means to perform public recitations of the Work and to communicate to the public those public recitations, by any means or process, including by wire or wireless means or public digital performances; to make available to the public Works in such a way that members of the public may access these Works from a place and at a place individually chosen by them; to perform the Work to the public by any means or process and the communication to the public of the performances of the Work, including by public digital performance; to broadcast and rebroadcast the Work by any means including signs, sounds or images.
 - i. **"Reproduce"** means to make copies of the Work by any means including without limitation by sound or visual recordings and the right of fixation and reproducing fixations of the Work, including storage of a protected performance or phonogram in digital form or other electronic medium.

2. Fair Dealing Rights. Nothing in this License is intended to reduce, limit, or restrict any uses free from copyright or rights arising from limitations or exceptions that are provided for in connection with the copyright protection under copyright law or other applicable laws.

3. License Grant. Subject to the terms and conditions of this License, Licensor hereby grants You a worldwide, royalty-free, non-exclusive, perpetual (for the duration of the applicable copyright) license to exercise the rights in the Work as stated below:

- a. to Reproduce the Work, to incorporate the Work into one or more Collections, and to Reproduce the Work as incorporated in the Collections;
- b. to create and Reproduce Adaptations provided that any such Adaptation, including any translation in any medium, takes reasonable steps to clearly label, demarcate or otherwise identify that changes were made to the original Work. For example, a

- translation could be marked "The original work was translated from English to Spanish," or a modification could indicate "The original work has been modified.;"
- c. to Distribute and Publicly Perform the Work including as incorporated in Collections; and,
 - d. to Distribute and Publicly Perform Adaptations.
 - e. For the avoidance of doubt:
 - i. **Non-waivable Compulsory License Schemes.** In those jurisdictions in which the right to collect royalties through any statutory or compulsory licensing scheme cannot be waived, the Licensor reserves the exclusive right to collect such royalties for any exercise by You of the rights granted under this License;
 - ii. **Waivable Compulsory License Schemes.** In those jurisdictions in which the right to collect royalties through any statutory or compulsory licensing scheme can be waived, the Licensor waives the exclusive right to collect such royalties for any exercise by You of the rights granted under this License; and,
 - iii. **Voluntary License Schemes.** The Licensor waives the right to collect royalties, whether individually or, in the event that the Licensor is a member of a collecting society that administers voluntary licensing schemes, via that society, from any exercise by You of the rights granted under this License.

The above rights may be exercised in all media and formats whether now known or hereafter devised. The above rights include the right to make such modifications as are technically necessary to exercise the rights in other media and formats. Subject to Section 8(f), all rights not expressly granted by Licensor are hereby reserved.

4. Restrictions. The license granted in Section 3 above is expressly made subject to and limited by the following restrictions:

- a. You may Distribute or Publicly Perform the Work only under the terms of this License. You must include a copy of, or the Uniform Resource Identifier (URI) for, this License with every copy of the Work You Distribute or Publicly Perform. You may not offer or impose any terms on the Work that restrict the terms of this License or the ability of the recipient of the Work to exercise the rights granted to that recipient under the terms of the License. You may not sublicense the Work. You must keep intact all notices that refer to this License and to the disclaimer of warranties with every copy of the Work You Distribute or Publicly Perform. When You Distribute or Publicly Perform the Work, You may not impose any effective technological measures on the Work that restrict the ability of a recipient of the Work from You to exercise the rights granted to that recipient under the terms of the License. This Section 4(a) applies to the Work as incorporated in a Collection, but this does not require the Collection apart from the Work itself to be made subject to the terms of this License. If You create a Collection, upon notice from any Licensor You must, to the extent practicable, remove from the Collection any credit as required by Section 4(b), as requested. If You create an Adaptation, upon notice from any Licensor You must, to the extent practicable, remove from the Adaptation any credit as required by Section 4(b), as requested.
- b. If You Distribute, or Publicly Perform the Work or any Adaptations or Collections, You must, unless a request has been made pursuant to Section 4(a), keep intact all copyright notices for the Work and provide, reasonable to the medium or means You are utilizing: (i) the name of the Original Author (or pseudonym, if applicable) if supplied, and/or if the Original Author and/or Licensor designate another party or parties (e.g., a sponsor institute, publishing entity, journal) for attribution ("Attribution Parties") in Licensor's copyright notice, terms of service or by other reasonable means, the name of such party or parties; (ii) the title of the Work if supplied; (iii) to the extent reasonably practicable, the URI, if any, that Licensor

- specifies to be associated with the Work, unless such URI does not refer to the copyright notice or licensing information for the Work; and (iv) , consistent with Section 3(b), in the case of an Adaptation, a credit identifying the use of the Work in the Adaptation (e.g., "French translation of the Work by Original Author," or "Screenplay based on original Work by Original Author"). The credit required by this Section 4 (b) may be implemented in any reasonable manner; provided, however, that in the case of a Adaptation or Collection, at a minimum such credit will appear, if a credit for all contributing authors of the Adaptation or Collection appears, then as part of these credits and in a manner at least as prominent as the credits for the other contributing authors. For the avoidance of doubt, You may only use the credit required by this Section for the purpose of attribution in the manner set out above and, by exercising Your rights under this License, You may not implicitly or explicitly assert or imply any connection with, sponsorship or endorsement by the Original Author, Licensor and/or Attribution Parties, as appropriate, of You or Your use of the Work, without the separate, express prior written permission of the Original Author, Licensor and/or Attribution Parties.
- c. Except as otherwise agreed in writing by the Licensor or as may be otherwise permitted by applicable law, if You Reproduce, Distribute or Publicly Perform the Work either by itself or as part of any Adaptations or Collections, You must not distort, mutilate, modify or take other derogatory action in relation to the Work which would be prejudicial to the Original Author's honor or reputation. Licensor agrees that in those jurisdictions (e.g. Japan), in which any exercise of the right granted in Section 3(b) of this License (the right to make Adaptations) would be deemed to be a distortion, mutilation, modification or other derogatory action prejudicial to the Original Author's honor and reputation, the Licensor will waive or not assert, as appropriate, this Section, to the fullest extent permitted by the applicable national law, to enable You to reasonably exercise Your right under Section 3(b) of this License (right to make Adaptations) but not otherwise.

5. Representations, Warranties and Disclaimer

UNLESS OTHERWISE MUTUALLY AGREED TO BY THE PARTIES IN WRITING, LICENSOR OFFERS THE WORK "AS IS" AND MAKES NO REPRESENTATIONS OR WARRANTIES OF ANY KIND CONCERNING THE WORK, EXPRESS, IMPLIED, STATUTORY OR OTHERWISE, INCLUDING, WITHOUT LIMITATION, WARRANTIES OF TITLE, MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, NON-INFRINGEMENT, OR THE ABSENCE OF LATENT OR OTHER DEFECTS, ACCURACY, OR THE PRESENCE OF ABSENCE OF ERRORS, WHETHER OR NOT DISCOVERABLE. SOME JURISDICTIONS DO NOT ALLOW THE EXCLUSION OF IMPLIED WARRANTIES, SO SUCH EXCLUSION MAY NOT APPLY TO YOU.

6. Limitation on Liability. EXCEPT TO THE EXTENT REQUIRED BY APPLICABLE LAW, IN NO EVENT WILL LICENSOR BE LIABLE TO YOU ON ANY LEGAL THEORY FOR ANY SPECIAL, INCIDENTAL, CONSEQUENTIAL, PUNITIVE OR EXEMPLARY DAMAGES ARISING OUT OF THIS LICENSE OR THE USE OF THE WORK, EVEN IF LICENSOR HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

7. Termination

- a. This License and the rights granted hereunder will terminate automatically upon any breach by You of the terms of this License. Individuals or entities who have received Adaptations or Collections from You under this License, however, will not have their licenses terminated provided such individuals or entities remain in full compliance with those licenses. Sections 1, 2, 5, 6, 7, and 8 will survive any termination of this License.
- b. Subject to the above terms and conditions, the license granted here is perpetual (for the duration of the applicable copyright in the Work). Notwithstanding the above, Licensor reserves the right to release the Work under different license terms or to stop distributing the Work at any time; provided, however that any such election will not serve to withdraw this License (or any

other license that has been, or is required to be, granted under the terms of this License), and this License will continue in full force and effect unless terminated as stated above.

8. Miscellaneous

- a. Each time You Distribute or Publicly Perform the Work or a Collection, the Licensor offers to the recipient a license to the Work on the same terms and conditions as the license granted to You under this License.
- b. Each time You Distribute or Publicly Perform an Adaptation, Licensor offers to the recipient a license to the original Work on the same terms and conditions as the license granted to You under this License.
- c. If any provision of this License is invalid or unenforceable under applicable law, it shall not affect the validity or enforceability of the remainder of the terms of this License, and without further action by the parties to this agreement, such provision shall be reformed to the minimum extent necessary to make such provision valid and enforceable.
- d. No term or provision of this License shall be deemed waived and no breach consented to unless such waiver or consent shall be in writing and signed by the party to be charged with such waiver or consent.
- e. This License constitutes the entire agreement between the parties with respect to the Work licensed here. There are no understandings, agreements or representations with respect to the Work not specified here. Licensor shall not be bound by any additional provisions that may appear in any communication from You. This License may not be modified without the mutual written agreement of the Licensor and You.
- f. The rights granted under, and the subject matter referenced, in this License were drafted utilizing the terminology of the Berne Convention for the Protection of Literary and Artistic Works (as amended on September 28, 1979), the Rome Convention of 1961, the WIPO Copyright Treaty of 1996, the WIPO Performances and Phonograms Treaty of 1996 and the Universal Copyright Convention (as revised on July 24, 1971). These rights and subject matter take effect in the relevant jurisdiction in which the License terms are sought to be enforced according to the corresponding provisions of the implementation of those treaty provisions in the applicable national law. If the standard suite of rights granted under applicable copyright law includes additional rights not granted under this License, such additional rights are deemed to be included in the License; this License is not intended to restrict the license of any rights under applicable law.

20. This software uses the Newtonsoft.Json 13.0.3 library

(<https://www.newtonsoft.com/json>)

This software makes use of the Newtonsoft.Json library. This library has not been modified by Akoya and all rights are reserved by the copyright holder. This software is provided "as is" without express or implied warranty, and with no claim as to its suitability for any purpose. The library is covered by the MIT License which is reproduced below.

Copyright © 2007 James Newton-King

The MIT License (MIT)

Permission is hereby granted, free of charge, to any person obtaining a copy of this software and associated documentation files (the "Software"), to deal in the Software without restriction, including without limitation the rights to use, copy, modify, merge, publish, distribute, sublicense, and/or sell copies of the Software, and to permit persons to whom the Software is furnished to do so, subject to the following conditions:

The above copyright notice and this permission notice shall be included in all copies or substantial portions of the Software.

THE SOFTWARE IS PROVIDED "AS IS", WITHOUT WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND NONINFRINGEMENT. IN NO EVENT SHALL THE AUTHORS OR COPYRIGHT HOLDERS BE LIABLE FOR ANY CLAIM, DAMAGES OR OTHER LIABILITY, WHETHER IN AN ACTION OF CONTRACT, TORT OR OTHERWISE, ARISING FROM, OUT OF OR IN CONNECTION WITH THE SOFTWARE OR THE USE OR OTHER DEALINGS IN THE SOFTWARE.

21. This software uses the CommandLine 1.6.0.0 library

(<https://github.com/commandlineparser/commandline>)

This software makes use of the Autofac library. This library has not been modified by Akoya and all rights are reserved by the copyright holder. This software is provided "as is" without express or implied warranty, and with no claim as to its suitability for any purpose. The library is covered by the MIT License.

Copyright © 2005 - 2012 Giacomo Stelluti Scala

22. This software uses the SendGrid 9.28.1 library (<https://sendgrid.com>)

This software makes use of the SendGrid library. This library has not been modified by Akoya and all rights are reserved by the copyright holder. This software is provided "as is" without express or implied warranty, and with no claim as to its suitability for any purpose. The library is covered by the MIT License.

Copyright © 2023, Twilio SendGrid, Inc. help@twilio.com

Should you have any question concerning this Agreement, you may contact Akoya by writing to Akoya Biosciences, 100 Campus Drive, 6th Floor, Marlborough, MA 01752. You may also call +1 855.896.8401

Phenomager HT

USER GUIDE

To learn more visit [AKOYABIO.COM](https://www.akoyabio.com) or email us at [INFO@AKOYABIO.COM](mailto:info@akoyabio.com)

For Research Use Only. Not for use in diagnostic procedures.

© 2024 Akoya Biosciences, Inc. All rights reserved. All trademarks are the property of Akoya Biosciences unless otherwise specified.

PD-000010 Rev B